



NEUROSURGERY

AESCULAP® NEUROENDOSCOPY

Intraventricular, Endoscope-Assisted,
Transnasal Neuroendoscopic Equipment

With comments from international experts in the field
of neuroendoscopy and minimally-invasive neurosurgery.


AESCULAP® NEUROENDOSCOPY

In 1924, the famous general and neurological surgeon William Halsted expressed his belief "...that the tendency will always be in the direction of exercising greater care and refinement in operating". Today, within the third millennium this fundamental philosophy of minimally invasive therapy should be emphasized more than ever before, operating with a minimum of iatrogenic trauma while achieving maximum surgical efficiency.

Recent improvements in preoperative imaging and surgical instrumentation allow neurosurgeons to treat more complex pathologies through customized less invasive approaches.

Using the advanced diagnostic tools of digital subtraction angiography, 3D angiography, computer tomography and magnetic resonance imaging, one

is able to demonstrate and elucidate preoperatively the individual anatomy and pathology of the patient. Therefore, anatomically preformed surgical dissection can be described preoperatively and may so be included into the planning of surgery. With the individual anatomic details of a specific patient, it becomes possible to perform a tailored surgical procedure reducing the size of the





Michael Fritsch
Neubrandenburg, Germany



Jeremy Greenlee
Iowa City, USA



André Grotenhuis
Nijmegen, Netherlands



Peter Nakaji
Phoenix, USA



Mark Souweidane
New York, USA



Charles Teo
Sydney, Australia

skin incision, the craniotomy, and the extent of brain surface traumatization and retraction to a necessary minimum limit. These advantages of minimally invasive microsurgery contribute to improved postoperative results, including shorter hospitalization time because of reduction of the risk for complications.


However, small sized minimally invasive approaches cause two important limitations: the significant loss of optical control and limited maneuverability of microsurgical instruments. The intraoperative use of endoscopes and dedicated minimally invasive instruments overcome these restrictions, thus enabling neurosurgeons to achieve deep seated regions without approach related traumatization of sensitive neurovascular structures.



The endoscopic image allows illumination and inspection of angles in hidden parts of the surgical field with the and clear depiction of anatomical details. In addition, due to the enormous optical depth of field of modern endoscopes, endoscopes provide a three dimensional aspect of anatomic structures. Recently, the intraoperative use of full high definition (HD) image quality offers a new area in endoscopic neurosurgery with an increased range of indications in minimally invasive neurosurgery.

There are three main indications of endoscopic neurosurgery: the intraventricular, transcranial and transnasal application. In this brochure, contemporary endoscopic equipment and instrumentation is presented in a comprehensive way. International experts in the field of minimally invasive and endoscopic neurosurgery comment the different applications, giving remarks with important tips and ideas, thus providing valuable instructions for the use of endoscopes in the field of minimally invasive neurosurgery.

Michael Fritsch, Neubrandenburg, Germany
Jeremy Greenlee, Iowa City, USA
André Grotenhuis, Nijmegen, Netherlands
Peter Nakaji, Phoenix, USA
Mark Souweidane, New York, USA
Charles Teo, Sydney, Australia







INTRAVENTRICULAR NEUROENDOSCOPY

MINOP[®]

Intraventricular Neuroendoscopic System





The genesis of endoscopic surgery within the ventricular compartment can be attributed to the development of small caliber rod lens optics, fiberoptic light transmission and dedicated instrumentation. Since the advent of intraventricular endoscopic surgery, neurosurgeons have applied the technology to treat a number of disorders. While the enthusiasm has been great and the full potential not yet realized, a major benefit to the patient has been proven for selected conditions. Most notably the treatment of non-communicating hydrocephalus, management of patients with pineal region tumors, fenestration of intracranial cysts, and removal of colloid cysts have all been shown to provide significant benefit and reduced morbidity compared with conventional treatment strategies.

The benefit in minimally invasive endoscopic procedures is analogous to that of any endoscopic procedure, namely minimal tissue disruption, enhanced visualization, improved cosmetic results, shorter hospital stay, and less surgical morbidity. The surgeon willing to utilize intraventricular endoscopic surgery is first responsible for attaining a considerable degree of familiarity with the technology, relevant anatomy, and the surgical procedures. Given the relative nascence of the field, the discipline is only now being commonly implemented in training programs. Hence, for those that have not had the opportunity to have endoscopic surgery as part of their formal training, it is strongly recommended that the surgeon participates in established practical courses in endoscopic neurosurgery, such as the courses from the Aesculap Academy.

Once fluent with the endoscopic equipment, more advanced procedures can be performed with greater familiarity and experience. It is anticipated with future generations of neurosurgeons that the endoscope will be an indispensable part of the neurosurgeon's armamentarium given the unmatched image resolution and minimally invasive qualities.

This foreseeable integration will expectantly be paralleled with continued evolution in compatible equipment to suit the needs of an expanding repertoire.

Few neurosurgical procedures demand a degree of familiarity with equipment as do neuroendoscopic techniques. This feature is somewhat explained by the recent introduction of the neuroendoscope as well as the delicate nature of the equipment. The basic components of any neuroendoscopic procedure include the endoscope and trocar, a camera with light source and monitor, as well as compatible instrumentation.

Charles Teo
Mark Souweidane



Charles Teo
Sydney, Australia




Mark Souweidane
New York, USA

MINOP[®]

Intraventricular Neuroendoscopic System - MINOP[®] Trocars

- Smooth tip of trocar for atraumatic insertion into the brain
- Single obturator for working channel enables insertion of the trocar, under visual control, with the endoscope
- Large depth scale on the outer shaft of the trocar
- Conical entry of working channel for intuitive insertion of instruments into trocar
- Attachment on top of trocar for improved handling and universal connection of peripheral devices



FF399R

MINOP[®] trocar,

Outer diameter 6 mm

4 channels:

- Endoscope channel: diam. 2.8 mm
- Working channel: diam. 2.2 mm
- Irrigation channel: diam. 1.4 mm
- Overflow channel, diam. 1.4 mm

Including 4 obturators for all channels



"I had used the Aesculap MINOP[®] system for all intraventricular cases and was mostly pleased with its versatility and safety. However, I had some concerns regarding its user-friendliness and applicability when one needed to be a 2-handed surgeon. Both these issues have been addressed with the new, improved MINOP[®] trocar and I have been very pleased with its added safety and practicality. I honestly believe it is quite clearly the best scope on the market for intraventricular endoscopic procedures. I applaud Aesculap for listening to the people who count most... the surgeons!"

Charles Teo, Sydney, Australia





FF398R

MINOP® trocar,



Outer diam. 4.6 mm

3 channels:

- Scope channel, diam. 2.8 mm
- Irrigation channel, diam. 0.8 mm
- Overflow channel, diam. 0.8 mm

Including one obturator for scope channel

1 sealing cap for pressure balance in endoscope channel



FF397R


MINOP® trocar,

outer diam. 3.2 mm

1 channel:

- Single channel for endoscope, including 1 obturator
- Endoscope channel: diam. 2.8 mm


1 sealing cap for pressure balance in scope channel



MINOP[®]

Intraventricular Neuroendoscopic System - MINOP[®] Endoscopes

- FULL HD compatible
- New optical components for enlarged image area and enhanced image quality, brightness, contrast
- Improved fibre optics provide more light
- The external tube is made from a high strength special alloy for superior breaking resistance
- Service-optimised construction reduces maintenance costs
- Autoclavable/Steris/Sterrad



180 mm, 7"

PE184A


MINOP[®] angled endoscope

Direction of view: 0° (green ring)

Shaft diam.: 2.7 mm

Shaft length, 180 mm

Autoclavable



180 mm, 7"

PE204A

MINOP[®] angled endoscope

Direction of view: 30° upwards (red ring)

Shaft diam.: 2.7 mm

Shaft length: 180 mm

Autoclavable



Instruments


- Shaft length 265 mm
- Diameter 2.0 mm
- Fully detachable for reprocessing
- High precision instrument tip

Tactile Feedback

- Integrated tactile feedback delivers small resistance indicating that instrument tip emerges from the trocar
- Improves control during insertion of instruments

Rotating Knob

- By rotating the knob slightly with index finger, the tip of instrument turns equally
- No need anymore to turn/rotate instrument with the entire arm/handle
- Improves precision of neuroendoscopic surgery
- Integrated safety mechanism in instrument shaft




"A very appealing feature of the MINOP tube shaft instruments is a rotational capability of the instrument tip through a coaxial system thus eliminating the need for hand rotation and reducing excessive movement of the endoscope. Irrespective of the instrument, graduated markings or precalibrated indicators on the shaft are important in providing the surgeon knowledge as to when the instrument will enter the endoscopic field. Even more safety is provided by the new tactile feedback of the improved MINOP instruments. A small spring delivers a tactile resistance "telling" the surgeon that the instrument tip is exiting the trocar."


Mark Souweidane, New York, USA

MINOP[®]

Intraventricular Neuroendoscopic System - MINOP[®] Rigid Instruments



Instrument complete: Handle · outer tube · jaw part with inner tube



FF385R

MINOP[®] scissors

sharp/sharp



FF387R

MINOP[®] biopsy forceps



FF386R

MINOP[®] scissors

blunt/blunt



FF388R

MINOP[®] grasping and dissecting forceps



FF389R

MINOP[®] surgical forceps, 1 x 2 teeth

The very delicate MINOP[®] instruments should be carefully detached completely and be pre-cleaned manually at the end of the operation. Keeping them in dedicated trays for reprocessing and sterilization protects the super-fine instrument tips. A careful handling by trained operating & CSSD staff is highly recommended and can eliminate the wear and tear of these sensitive but highly necessary neuroendoscopic tools.

MINOP[®]

Intraventricular Neuroendoscopic System - MINOP[®] Rigid Instruments - Spare Parts

 2 mm Jaw part with inner tube for FF385R - FF389R

FF433R

MINOP[®] outer tube, only

FF432R

MINOP[®] instrument handle, only



FF435R

MINOP[®] scissors, jaw part

sharp/sharp



FF437R

MINOP[®] biopsy forceps, jaw part



FF436R

MINOP[®] scissors, jaw part

blunt/blunt




FF438R

MINOP[®] grasping and dissecting forceps, jaw part



FF439R

MINOP[®] surgical forceps, jaw part,
1 x 2 teeth



■ For disassembly and assembly of MINOP[®] tube shaft instruments see brochure no. C60911

MINOP®

Intraventricular Neuroendoscopic System - MINOP® Flexible Instruments

Ø 1 mm 1.0 mm Instruments, non detachable

250 mm, 10"



FF373R
MINOP® micro scissors



FF374R
MINOP® micro grasping and dissecting forceps



FF378R
MINOP® micro biopsy forceps

FF373R - FF378

With irrigation port for reprocessing/cleaning

Ø 1 mm 1.0 mm instruments, for bi-instrumental work

Flexible instruments:

- For bi-instrumental/bi-manual neuroendoscopic surgery
- E.g. grasping and cutting, grasping and coagulating, grasping and fenestrating
- To be used through irrigation or overflow channel of the MINOP® trocar FF399R
- Diam. 1.0 mm, shaft length 250 mm
- Non-detachable
- With irrigation port for reprocessing/cleaning



"The MINOP® system is providing bi-instrumental endoscopic work. For example in cyst removal or endoscopic tumor surgery the surgeon has the opportunity to grasp and cut or grasp and coagulate at the same time. One can utilize flexible instruments or electrodes in one of the side-channels and rigid tube shaft instruments in the working channel. The design of the side-channels of the MINOP® trocar makes sure that both instruments do not interfere with each other."

Michael Fritsch, Neubrandenburg, Germany



MINOP[®]

Intraventricular Neuroendoscopic System - MINOP[®] Electrodes

Monopolar electrodes

for bi-instrumental / bi-manual working through irrigation/overflow channel of trocar FF399R

GK361R

Blunt electrode, diam. 1.1 mm



GK363R

Needle electrode, diam. 1.1 mm



Not for bi-manual

GK364R

Hook electrode, 45°, diam. 2.2 mm



GK365R

Hook electrode, 70°, diam. 2.2 mm



GK362R

Hook electrode, 90°, diam. 2.2 mm



GK366R

Hook electrode, 180°, diam. 2.2 mm



GN202

Monopolar cable, 3.5 m length
suitable for GN300, GN640



Bipolar electrode

GK360R

Fork electrode, 0°, diam. 2.1 mm



GN130

Bipolar cable, 4 m length
suitable for GN060, GN160, GN300,
GN640



- For further details see brochures no. C46702 and no. C61511.



MINOP[®]

Intraventricular Neuroendoscopic System - MINOP[®] Single use Suction Cannulas

- For removal of cystic intraventricular lesions
- For puncturing the floor of the 3rd ventricle
- Depth scale
- Outer diameter of 2.0 mm
- Suitable for working channel of MINOP[®] trocar FF399R
- Available with blunt or sharp tip suction cannula
- Optional control of suction
 - via thumb plate or
 - via syringe
- Single-use, sterile packed



FH606SU

Suction cannula

blunt tip 0°, diam. 2 mm



FH607SU

Suction cannula


sharp tip 45°, diam. 2 mm



MINOP[®]

Intraventricular Neuroendoscopic System - MINOP[®] Single use Introducer

- Introducer sheath protects the brain while inserting and removing the endoscope/trocar
- Especially for MINOP[®] trocar FF399R
- 19 Fr disposable introducer set including obturator and sheath
- Round & blunt obturator tip for atraumatic insertion into the ventricles
- Large depth scale
- Easy to peel with side handles
- Sterile packed



FH604SU

Introducer

19 Fr

Sales unit:


PAK - Package of 5 pieces

The MINOP[®] suction cannula and the MINOP[®] disposable introducer can be used in almost any intraventricular neuroendoscopic surgery providing more control during the procedure. The suction cannula can be used for the controlled and fast removal of intraventricular soft tumors or colloid cysts with its sharp cannula tip or even for the opening of the floor of the 3rd ventricle. The disposable introducer (also called peel away) is very helpful when several intraparenchymal in- and out-movements of the trocar are necessary.

MINOP®

Intraventricular Neuroendoscopic System - MINOP® Storage

■ for MINOP® trocars and endoscopes




FF358R Dimension (L/W/H) 489x257x63 mm

Storage rack with silicone protection cushioning tray and lid

only for reprocessing, not for transportation/shipment
(instruments not included)

■ for MINOP® instruments and electrodes



FF359R Dimension (L/W/H) 485x253x120 mm

Storage rack with silicone protection cushioning tray without lid (lid not necessary)

only for reprocessing, not for transportation/shipment
(instruments not included)

■ 1/1 Container (basic version) for storage racks FF358R and FF359R



consisting of:

JK440

Container body 1/1 for FF358R
without base perforation

Outside/Inside dimensions with lid:

L/W/H 592 x 285 x 108 mm

L/W/H 544 x 258 x 75 mm

JK444

Container body 1/1 for FF359R
without base perforation

Outside/Inside dimensions with lid:

L/W/H 592 x 285 x 209 mm

L/W/H 544 x 258 x 172 mm

JK486

Inner lid 1/1

blue

Dedicated storage racks for cleaning and reprocessing are highly recommended for your neuroendoscopic equipment. A special-designed storage concept is keeping the scopes and instruments properly stored and protected.



PF893800

Cleaning brush



EJ751251

Sealing cap Luer-Lock

Sales unit:

PAK = Package of 20 piece

EJ751200

Sealing cap Luer-Lock

Sales unit:

ST = Package of 1 piece



■ For further details see brochure no. C40402

PAEDISCOPE®

Paediatric Intraventricular Neuroendoscopic System

PAEDISCOPE®

- 3.0 mm outer diameter for minimally invasive paediatric surgery
- 30.000 pixel fibre optic
- Light-weight and ergonomic design
- Weight of camera is away from the operating site
- Depth scale



150 mm, 6"

PF010A

Endoscope shaft
with integrated optical fibres



PF011A


Ocular with focus
for complete PAEDISCOPE® please order
both: PF010A and PF011A



PAEDISCOPE®

Paediatric Intraventricular Neuroendoscopic System - Flexible Instruments

 1 mm Instruments, non detachable




FF373R
Micro scissors

FF374R
Micro grasping and dissecting forceps

FF378R
Micro biopsy forceps

FF373R - FF378
With irrigation port for reprocessing/cleaning



FH603SU

PAEDISCOPE® Single use introducer, 10 FR

Sales unit:

PAK = Package of 5 pieces

- Introducer sheath protects the brain while inserting and removing the endoscope/trocar
- Especially for PAEDISCOPE® PF010A
- 10 Fr disposable introducer set including obturator and sheath
- Round & blunt obturator tip for atraumatic insertion into the ventricles
- Large depth scale
- Easy to peel with side handles
- Sterile packed



The peel away sheath protects the brain while inserting and removing the pediatric endoscope. Because of its small outer diameter, the Paediscop does not have a dedicated trocar. The blunt obturator tip of the sheath allows atraumatic insertion into the ventricles. The sheath has a depth scale for precise positioning and is easy to peel back the side handles. Using a peel away sheath is especially helpful, if repeated in and out movements of the scope are necessary or different instruments or catheters (e.g. for aqueductoplasty) have to be utilized in addition to the scope.

Michael Fritsch, Neubrandenburg, Germany


PAEDISCOPE®

Paediatric Intraventricular Neuroendoscopic System - Monopolar Electrodes



GK361R

Blunt electrode, diam. 1.1 mm



GK363R

Needle electrode, diam. 1.1 mm



GK202

Monopolar cable, length 3.5 m
suitable for GN300, GN640

PAEDISCOPE®

Paediatric Intraventricular Neuroendoscopic System - Storage

- Storage rack for PAEDISCOPE® instruments and electrodes



FF379R Dimension (L/W/H) 489 x 257 x 63 mm

Storage rack with silicone protection cushioning, tray and lid

only for reprocessing, not for transportation/shipment (instruments not included)

- 1/1 Container (basic version) for storage rack FF379R



consisting of:

JK440

Container body 1/1 for FF358R without base perforation

Outside/Inside dimensions with lid:

L/W/H 592 x 285 x 108 mm

L/W/H 544 x 258 x 75 mm

JK486


Inner lid 1/1

blue



- For further details see brochure no. C29202 and no. C40402







MINOP[®] InVent

Intraventricular Neuroendoscopic System for Experts - Trocar

- Larger trocar with an oval working channel allows innovative treatment options and multi-directional flexibility
- Bi-instrumental technique is similar to traditional micro neurosurgery, due to the increased freedom of movement
- For the first time, angled instruments can be used
- Up to 32 different instruments usable



FH620R

MINOP[®] InVent trocar,

Outer diameter 8,3 mm

3 (4) channels:

- Endoscope channel: diam. 2.8 mm
- Irrigation channel: diam. 2.2 mm

2 merging channels

- Large working/overflow channel: 3.7 x 6.5 mm
- Small working channel, diam. 2.2 mm


including 2 obturators for endoscope channel and large working/overflow channel



MINOP[®] InVent

Intraventricular Neuroendoscopic System for Experts - Endoscope

- FULL HD compatible
- New optical components for enlarged image area and enhanced image quality, brightness, contrast
- Improved fibre optics provide more light
- The external tube is made from a high strength special alloy for superior breaking resistance
- Service-optimised construction reduces maintenance costs
- Autoclavable/Steris/Sterrad



PE204A

MINOP[®] InVent angled endoscope

Direction of view: 30° upwards (red ring)

Shaft diam.: 2.7 mm

Autoclavable




"The MINOP[®] InVent system is truly unique and the next step for the future of Neuroendoscopy. This system allows for a true bi-manual technique through the large/small working channels expanding the possibilities to treat further indications. The angled instrumentation provide the ability to simultaneously grasp and cut or grasp and coagulate similar to traditional microsurgery. The MINOP InVent provides a new possibility for the treatment of intra- and paraventricular cysts and tumors in complex hydrocephalus and alleviating the need for certain craniotomies."

Mark Souweidane, New York, USA

MINOP[®] InVent

Intraventricular Neuroendoscopic System for Experts - Dissectors, Hook and Knife



FH629R

MINOP[®] InVent dissector,
tip width 2.2 mm



FH623R

MINOP[®] InVent hook, 90° blunt,
hook deflection width 3.5 mm



FH630R

MINOP[®] InVent dissector,
tip width 1.7 mm



FH623R

MINOP[®] InVent knife, backwards cutting,
knife deflection width 3.0 mm



FH631R

MINOP[®] InVent dissector,
tip width 1.0 mm





MINOP[®] InVent

Intraventricular Neuroendoscopic System for Experts - Shaft Instruments



Width x Height:
2.0 mm x 3.1 mm



FH621R

MINOP[®] InVent forceps, straight

FH622R

MINOP[®] InVent scissors, straight



FH622R

MINOP[®] InVent forceps, right



FH626R

MINOP[®] InVent scissors, left



FH623R

MINOP[®] InVent forceps, left



FH627R

MINOP[®] InVent scissors, right



FH624R

MINOP[®] InVent grasping forceps, straight



FH628R

MINOP[®] InVent scissors, upwards




MINOP[®] InVent


Intraventricular Neuroendoscopic System for Experts - Tube Shaft Instruments



2.0 mm **Instrument complete:** Handle · outer tube · jaw part with inner tube



- High precision working ends
- Rotation wheel for comfortable rotation of working end



FH635R

MINOP[®] InVent scissors
sharp/sharp



FH636R

MINOP[®] InVent scissors
blunt/blunt



FH636R

MINOP[®] InVent biopsy forceps



FH638R

MINOP[®] InVent grasping and dissecting forceps




FH639R

MINOP[®] InVent surgical forceps
1 x 2 teeth



MINOP[®] InVent

Intraventricular Neuroendoscopic System for Experts – Tube Shaft Instruments – Spare Parts



FH635200

MINOP[®] InVent outer tube, only



FH435R

MINOP[®] InVent scissors, jaw part
sharp/sharp



FH436R

MINOP[®] InVent scissors, jaw part
blunt/blunt



FH438R

MINOP[®] InVent grasping and
dissecting forceps, jaw part



FH633R

MINOP[®] InVent instrument
handle, only




FH436R

MINOP[®] InVent biopsy forceps, jaw part



FH439R

MINOP[®] InVent surgical forceps, jaw part
1 x 2 teeth





- For disassembly and assembly of MINOP[®] tube shaft instruments see brochure no. C60911

MINOP[®] InVent

Intraventricular Neuroendoscopic System for Experts - Bipolar Electrodes



Width x Height:
3.2 mm x 2.1 mm



GK343R

MINOP[®] InVent bipolar electrode,
0°, diam.: 2.7 mm



GK344R

MINOP[®] InVent bipolar electrode,
40°, diam.: 2.7 mm




GK345R

MINOP[®] InVent bipolar electrode,
30°, diam.: 2.7 mm



2.1 mm



GK360R

Fork electrode, 0°, diam. 2.1 mm


GN130

Bipolar cable, 4 m length
suitable for GN060, GN160, GN300,
GN640



MINOP[®] InVent

Intraventricular Neuroendoscopic System for Experts - Monopolar Electrodes



GK361R
Blunt electrode, diam. 1.1 mm



GK364R
Hook electrode, 45°, diam. 2.2 mm



GK363R
Needle electrode, diam. 1.1 mm



GK365R
Hook electrode, 70°, diam. 2.2 mm



GK362R
Hook electrode, 90°, diam. 2.2 mm






GK366R
Hook electrode, 180°, diam. 2.2 mm

GN202
Monopolar cable, 3.5 m length
suitable for GN300, GN640



MINOP[®] InVent

Intraventricular Neuroendoscopic System for Experts - Flexible Instruments,
Suction Cannulas



FH606SU
Suction cannula
blunt tip 0°, diam. 2 mm




FH607SU
Suction cannula
sharp tip 45°, diam. 2 mm



MINOP[®] InVent

Intraventricular Neuroendoscopic System for Experts - Storage

- for MINOP[®] InVent trocars and endoscopes



FH358R Dimension (L/W/H) 540 x 253 x 56 mm

Storage rack with silicone protection cushioning tray and lid only for reprocessing, not for transportation/shipment (instruments not included)

- 1/1 Container (basic version) for storage racks FF358R and FF359R,



consisting of:

JK440

Container body 1/1 for FF358R without base perforation

Outside/Inside dimensions with lid:

L/W/H 592 x 285 x 108 mm

L/W/H 544 x 258 x 75 mm

JK444


Container body 1/1 for FF359R without base perforation

Outside/Inside dimensions with lid:

L/W/H 592 x 285 x 209 mm

L/W/H 544 x 258 x 172 mm

- for MINOP[®] InVent instruments and electrodes



FH359R Dimension (L/W/H) 540 x 253 x 120 mm

Storage rack with silicone protection and cushioning, tray without lid (lid not necessary)

only for reprocessing, not for transportation/shipment (instruments not included)


JK486

Inner lid 1/1

blue







ENDOSCOPE-ASSISTED MICRO NEUROSURGERY



MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery



"The aim of minimally invasive neurosurgery is to avoid approach-related traumatization of the patient by creating a tailor-made limited craniotomy based on skilled preoperative planning.

Using modern diagnostic tools, surgical instruments and visual equipment, the specific anatomy and pathology of the individual patient can be precisely visualized and anatomical pathways and surgical corridors determined for the surgical approach. According to the predefined access, surgical dissection can be subsequently performed creating a much less traumatic cranial opening. The aim is not the limited cranial opening, but the limited approach associated injury with less brain exploration and retraction. The craniotomy should be as small as possible for minimally invasive exposure, but as large as necessary for achieving maximal surgical effect. In this way, limited exposure is not the primary goal but the result of the keyhole concept with the main and most important goal being to avoid surgery-related complications.

The intraoperative use of microscopes is mandatory in keyhole neurosurgery. The operating microscope provides both stereoscopic magnification and illumination of the surgical field. However, the loss of light intensity in the depth of the surgical field is a fundamental problem in keyhole approaches. For the purpose of bringing light into the site, operating microscopes can effectively be combined with the intraoperative use of modern endoscopes. The advantages of the endoscopic image are increased light, extended viewing angle and a better depiction of anatomical details in close-up. The endoscope

is especially ideal for obtaining a detailed view of structures in the shadow of the microscope's light beam. Thus, in situations during microsurgical dissection where additional visual information of the target area is desired or when avoidance of retraction of superficial structures is recommended, an endoscope may be introduced into the surgical site.

The use of dedicated microneurosurgical instruments is obligatory in transcranial endoscope-assisted microneurosurgery. Highly sophisticated instrumentation including microdrills, KERRISON micropunches, self-retaining retractors, suction tubes, fine bipolar forceps, microscissors, diamond knives, microforceps, microdissectors, micro-curettes, and clip applicators are mandatory for microsurgical dissection.

All before mentioned surgical tools - the microscope, endoscope and dedicated surgical instruments - complement each other and contribute in a TEAM-work manner to the goal of the keyhole concept: the achievement of the smallest iatrogenic trauma with the highest therapeutic effect for the patients."

Peter Nakaji



Peter Nakaji, Phoenix, USA

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery - Angled "PERNECZKY" Endoscopes

- FULL HD compatible
- Different viewing angles (0°, 30°, 70°)
- Angled endoscope design and lateral connection for camera and light source allow to use microscope and micro instruments in parallel
- Ergonomic handling by centered balance of weight
- Autoclavable/Steris[®]/Sterrads[®]
- Robust and rigid scope sheath




PE486A

Angled endoscope

Direction of view: 0°

Shaft diameter: 4 mm

Shaft length: 150 mm, 6"



"I have been using the Aesculap angled Perneczky scopes since the mid nineties and in over 1000 cases. I have trialed many different scopes for endoscope-assisted surgery but the Perneczky scopes have the versatility that I need when removing tumors from many different cranial locations. The main advantage of the angled scopes is the unique design that allows simultaneous use of endoscope and microscope. Other important qualities that are met by this system are robustness, ability to use it to retract if necessary and clarity of image. I believe these scopes are an essential tool in the neurosurgeon's armamentarium."

Charles Teo, Sydney, Australia




PE506A**Angled endoscope**

Direction of view: 30°, upwards

Shaft diameter: 4 mm


Shaft length: 150 mm, 6"

**PE526A****Angled endoscope**

Direction of view: 70°, upwards

Shaft diameter: 4 mm

Shaft length: 150 mm, 6"

**JF324R****Storage tray**

with silicone cushioning racks and lid
for 2 angled neuroscopes (not included)
(L/W/H 247 x 257 x 64 mm)



"During microneurosurgical skull base approaches for either vascular lesions or tumors, there is often a difficulty of visualizing important neurovascular structures around and behind the lesion. In such a situation, the use of endoscopes has greatly advanced my surgical possibilities. The additional view through the endoscopes, which is complementary to what can be seen through the operating microscope, facilitates the handling of the lesion, be it aneurysm clipping or tumor removal, while at the same time there is no need for extensive retraction or bone removal."

André Grotenhuis, Nijmegen, Netherlands


MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery – Aesculap Micro Instruments

Small craniotomies or narrow operative sites require especially designed fine and slender micro instruments

Experience our three different lines of minimally invasive Micro Instruments

page 43 –48



For more information about XS Micro Instruments please see our brochure C77011

page 49 – 58



For more information about MIN Set please see our brochure C92011



For more information about SENSATION Micro Instruments please see our brochure C84902

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery
XS Tube Shaft Micro Instruments

XS MICRO INSTRUMENTS


- Narrow tubular shaft for maximal freedom of movement
- Angled bayonet shape for enhanced sight lines and easier handling
- Dismountable for easy, effective cleaning and reprocessing
- Exchangeable handles and jaw inserts for economic use




MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery

XS Micro Scissors acc. PERNECZKY/CHRISTANTE, Bayonet-shaped



	Instrument complete	Single parts		Instrument complete	Single parts		Instrument complete	Single parts	
		Jaw insert	Handle		Jaw insert	Handle		Jaw insert	Handle
straight jaws									
sharp/sharp	FM670R	FM675R	FM730R	FM671R	FM676R	FM731R	FM672R	FM677R	FM732R
blunt/blunt	FM690R	FM695R		FM691R	FM696R		FM692R	FM697R	
Working length	70 mm, 2 3/4"			100 mm, 4"			130 mm, 5 1/8"		
Total length	200 mm, 8"			230 mm, 9"			260 mm, 10 1/4"		




	Instrument complete	Single parts		Instrument complete	Single parts		Instrument complete	Single parts	
		Jaw insert	Handle		Jaw insert	Handle		Jaw insert	Handle
curved jaws									
sharp/sharp	FM680R	FM685R	FM730R	FM681R	FM686R	FM731R	FM682R	FM687R	FM732R
blunt/blunt	FM700R	FM705R		FM701R	FM706R		FM702R	FM707R	
Working length	70 mm, 2 3/4"			100 mm, 4"			130 mm, 5 1/8"		
Total length	200 mm, 8"			230 mm, 9"			260 mm, 10 1/4"		





MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery -
XS Micro Tissue Forceps acc. PERNECZKY/CHRISTIANE, Bayonet-shaped



Jaw	Instrument complete	Single parts		Instrument complete	Single parts		Instrument complete	Single parts	
		Jaw insert	Handle		Jaw insert	Handle		Jaw insert	Handle
0.9 mm	FM710R	FM715R	FM730R	FM711R	FM716R	FM731R	FM712R	FM717R	FM732R
Working length	70 mm, 2 3/4"			100 mm, 4"			130 mm, 5 1/8"		
Total length	200 mm, 8"			230 mm, 9"			260 mm, 10 1/4"		



Jaw	Instrument complete	Single parts		Instrument complete	Single parts		Instrument complete	Single parts	
		Jaw insert	Handle		Jaw insert	Handle		Jaw insert	Handle
3 mm, sharp 	FM720R	FM725R	FM730R	FM721R	FM726R	FM731R	FM722R	FM727R	FM732R
Working length	70 mm, 2 3/4"			100 mm, 4"			130 mm, 5 1/8"		
Total length	200 mm, 8"			230 mm, 9"			260 mm, 10 1/4"		

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery XS Tube Shaft Aneurysm Clip Applying Forceps

- 360° rotation
- suitable for narrow approaches



	Standard	Mini
For Titanium clips	FT494T	FT489T
For Phynox clips	FE494K	FE489K



	Standard	Mini
For Titanium clips	FT495T	FT490T
For Phynox clips	FE495K	FE490K



	Standard	Mini
For Titanium clips	FT496T	FT491T
For Phynox clips	FE496K	FE491K

The cause for the significant superiority of the endovascular treatment of aneurysms compared with the surgical therapy in the ISAT study was the surgical morbidity and mortality of large sized standard approaches. In my opinion, surgical clipping will play an important role in the treatment of intracranial aneurysms in the future only, if it will be able to reduce approach related complications using limited craniotomies. The use of endoscope-assisted techniques and tube-shaft clip appliers offer increased safety in keyhole vascular neurosurgery, thus achieving the basic goal with minimally invasive and maximal effective aneurysm closure.

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery
MIN Noir[®] Micro Instruments



Slender design and angled bayonet shape

Improved visibility of the surgical site due to the slender design. Angled bayonet shape allows for less obstructions while working under the microscope.



Round golf ball handle design

Designed to provide an excellent grip and enable easy rotation of the instruments. This allows precise handling.



Various working lengths

One handle design aligned with precisely adapted working lengths. Always provides the right instrument at your fingertips!



Noir[®], No Irritating Reflections

Aesthetic surface coating effectively prevents disturbing light reflections.



Fine instrument tips

Especially important when working in very small operating corridors and close to sensitive structures.



MINOP[®] TEAM


Transcranial Endoscope Assisted Microneurosurgery

MIN Noir[®] Micro Scissors with Round Golfball Design Handle, Bayonet-shaped




	straight jaws	slightly curved jaws	curved jaws	straight jaws	slightly curved jaws	curved jaws	straight jaws	slightly curved jaws	curved jaws
--	---------------	----------------------	-------------	---------------	----------------------	-------------	---------------	----------------------	-------------

sharp/sharp



FD701B	FD702B	FD703B	FD731B	FD732B	FD733B	FD771B	FD771B	FD773B
--------	--------	--------	--------	--------	--------	--------	--------	--------

blunt/blunt



FD704B	FD705R	FD706B	FD734B	FD735B	FD736B	FD774B	FD775B	FD776B
--------	--------	--------	--------	--------	--------	--------	--------	--------

Working length

70 mm, 2 3/4"

100 mm, 4"

130 mm, 5 1/8"

Total length


200 mm, 8"

230 mm, 9"

260 mm, 10 1/4"




■ For more information to MIN Instruments, please ask your local Aesculap sales representative or see our brochure C92011



MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery

MIN Noir[®] Micro Needle Holders with Round Golfball Design Handle, Bayonet-shaped



	FD717B	FD718B	FD719B
Working length	70 mm, 2 3/4"	90 mm, 3 1/2"	120 mm, 4 3/4"
Total length	200 mm, 8 7/8"	220 mm, 8 3/4"	250 mm, 9 3/4"



MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery - MIN Tissue and Tumor Grasping Forceps



0.5 mm 	FD711B	FD741B	FD761B
0.9 mm 		FD743B	FD763B
Working length	70 mm, 2 3/4"	90 mm, 3 1/2"	120 mm, 4 3/4"
Total length	190 mm, 7 1/2"	210 mm, 8 1/4"	240 mm, 9 1/2"





2.5 mm 	FD766B	FD786B	FD767B	FD787B	FD768B	FD788B
3.5 mm 	FD769B	FD789B	-	-	-	-
Working length	90 mm, 3 1/2"	120 mm, 4 3/4"	90 mm, 3 1/2"	120 mm, 4 3/4"	90 mm, 3 1/2"	120 mm, 4 3/4"
Total length	210 mm, 8 1/4"	240 mm, 9 1/2"	210 mm, 8 1/4"	240 mm, 9 1/2"	210 mm, 8 1/4"	240 mm, 9 1/2"

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery
MIN Micro Instruments

		Noir [®] Modular Handles	
	FD811B	Handle, 8 mm, 100 mm	
	FD812B	Handle, 11 mm, 100 mm	
	FD818B	Handle, 8 mm, 100 mm	
	FD819B	Handle, 11 mm, 100 mm	
	FD848B	Handle, 8 mm, 100 mm	
	FD849B	Handle, 11 mm, 100 mm	

1/2

		Noir [®] Probes / Hooks	
	FD797B	Probe ball-tip, 200 mm, 0°	
	FD798B	Probe ball-tip, 200 mm, 45°	
	FD799B	Probe ball-tip, 200 mm, 90°	
	FD808B	Hook, blunt, 200 mm, 45°	
	FD809B	Hook, blunt, 200 mm, 90°	
	FD805B	Hook, sharp, 200 mm, 90°	

1/1

		Noir [®] Scoops	
	FD814B	Scoop, 200 mm, 2 mm, 10°	
	FD815B	Scoop, 200 mm, 2 mm, 45°	
	FD816B	Scoop with neck, 200 mm, 2 mm, 45°	

1/1

		Noir [®] Dissectors	
	FD821B	Dissector, curved, 200 mm, 1 mm	
	FD822B	Dissector, curved, 200 mm, 2 mm	
	FD823B	Dissector, curved, 200 mm, 3 mm	

1/1

Noir® Currettes

	FD824B	Curette, 200 mm, diam. 4 mm, 0°
	FD825B	Curette, 200 mm, diam. 4 mm, 45°
	FD826B	Curette, 200 mm, diam. 4 mm, 90°
	FD827B	Curette with neck, 200 mm, diam. 4 mm, 45°
	FD828B	Curette with neck, 200 mm, diam. 4 mm, 90°
	FD835B	Curette, 200 mm, diam. 6,5 mm, 45°
	FD836B	Curette, 200 mm, diam. 6,5 mm, 90°

1/1

Noir® Rasparatories


	FD831B	Rasparatory, 200 mm, 1 mm
	FD832B	Rasparatory, 200 mm, 2 mm
	FD833B	Rasparatory, 200 mm, 3 mm

1/1

Noir® Tumor knives

	FD839B	Noir Tumor knife, 200 mm, diam. 1,5 mm, 45°
	FD840B	Noir Tumor knife, 200 mm, diam. 3 mm, 45°
	FD841B	Noir Tumor knife, 200 mm, diam. 4,5 mm, 45°

1/1

**FD467R** Tray

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery
MIN Pivot-Point Bipolar Forceps



Working length



Aesculap tab connector 

0.7 mm  1/1

1.0 mm  1/1

	GK822R	GK826R	GK823R	GK827R
	GK824R	GK828R	GK825R	GK829R
Working length	95 mm, 3 3/4"	135 mm, 5 1/4"	95 mm, 3 3/4"	135 mm, 5 1/4"
Total length	215 mm, 8 1/2"	255 mm, 10"	215 mm, 8 1/2"	255 mm, 10"

"The black „pivot“ bipolar forceps are a great advance. The bipolar is as essential a tool as the neurosurgeon's own fingers. As we go more and more minimally invasive, the need for a very slim, responsive bipolar that will work under tight conditions is essential. The tips can be precisely separated even when the shafts are together in a tiny space. This is a must-have instrument, especially for transphenoidal and keyhole approaches."


Peter Nakaji, Phoenix, USA




MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery – Bipolar Yasargil Forceps

Bipolar Yasargil forceps:
extra-small bipolar forceps for
keyhole approaches



Aesculap tab connector 


	-	0.4 mm  1/1	GK780R	0.4 mm	-
0.7 mm  1/1	GK777R	0.7 mm  1/1	GK781R	0.7 mm  1/1	GK785R
Working length	95 mm, 3 3/4"		95 mm, 3 3/4"		95 mm, 3 3/4"
Total length	215 mm, 8 1/2"		215 mm, 8 1/2"		215 mm, 8 1/2"



- For more information please ask your local Aesculap sales representative or see our brochure C30481

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery
MIN Modular Suction Cannulas



GF005B

Thumb control for modular suction tubes, tear drop shaped for a very precise suction regulation, incl. mandrin.

	○S	○○M	○○○L	○○○○XL	
Working length	80 mm, 3 1/2"	100 mm, 4"	120 mm, 4 3/4"	140 mm, 5 1/2"	
4 Fr. straight	○	GF025B	GF035B	GF045B	GF055B
6 Fr. straight	○	GF026B	GF036B	GF046B	GF056B
6 Fr. straight, lateral holes	◊	-	GF038B	GF048B	GF058B
6 Fr. curved left	○	-	GF030B	-	-
6 Fr. curved right	○	-	GF031B	-	-
8 Fr. straight	◎	GF027B	GF037B	GF047B	GF057B
8 Fr. straight, lateral holes	◎	-	GF039B	GF049B	GF059B
8 Fr. curved left	◎	-	GF032B	-	-
8 Fr. curved left	◎	-	GF033B	-	-

3 Fr. = 1 mm


Benefits...

- Atraumatic tips enable blunt dissection and retraction
- Lateral bore holes reduce suction pressure
- Damaged suction tubes can be replaced individually
- Tray weight is reduced and less space is required
- Greater freedom for even more flexible use



MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery
RAABE Micro Suction Cannulas



Endoscope-Assisted

RAABE	○S	○○M	○○○L	○○○○XL	
Working length	80 mm, 3 1/8"	100 mm, 4"	120 mm, 4 3/4"	140 mm, 5 1/2"	
Total length	130 mm, 5 1/8"	150 mm, 6"	165 mm, 6 1/2"	185 mm, 7 1/4"	
4F yellow, 1.4 mm	○	GF470R	GF473R	GF476R	GF479R
6 F blue, 2.0 mm	◎	GF471R	GF474R	GF477R	GF480R
8F green, 2.7 mm	◎	GF472R	GF475R	GF478R	GF481R

3 Fr = 1 mm



The ball tip at the end of the instrument allows gentle preparation and stable atraumatic retraction.



Colour coding for rapid identification of all three diameters. Black Rings as indicators to identify the instrument length.

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery TREND Currettes and Dissectors

Bayonet design with ergonomic grasping part and semi-sharp tips for pituitary and skull base



FA041R-FA068R
Working length:
130 mm, 5 1/8"
Total length:
280 mm, 11"

diam. 6.5 mm		diam. 6.5 mm			
					
NICOLA FA041R	NICOLA FA042R	HARDY FA043R	HARDY FA044R		
Curette 45° vertical angled long neck	Curette 45° horizontal angled short neck	Enucleator left cutting	Enucleator right cutting		
diam. 4.0 mm		diam. 4.0 mm			
					
HARDY FA045R	HARDY FA046R	HARDY FA047R	HARDY FA060R		
Curette 90° left angled long neck	Curette 90° left angled short neck	Curette 90° right angled long neck	Curette 90° right angled short neck		

diam. 4.0 mm



HARDY
FA061R

Curette
45° horizontal,
left angled.
short neck

diam. 4.0 mm



HARDY
FA062R

Curette
45° horizontal,
right angled
short neck

diam. 6.0 mm



HARDY
FA063R

Curette
90° left angled
long neck

diam. 6.0 mm



HARDY
FA064R

Curette
90° left angled
short neck

diam. 6.0 mm



HARDY
FA065R

Curette
90° right angled
long neck

diam. 6.0 mm



HARDY
FA066R

Curette
90° right angled
short neck

diam. 1.7 mm



REULEN-
LANDOLT

FA067R
Micro Hook

diam. 2.0 mm



REULEN-
LANDOLT

FA068R

Dissector
blunt



MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery Diamond Knives

Diamond knives

Due to the properties of their unique blade material, diamond knives offer excellent cutting characteristics for highly clean, precise and force-free incisions especially in neurosurgical applications and in cardiovascular and thoracic surgery.



- Blade made of natural diamond
- Superior mechanical stability & elasticity of the blade
- Sustained sharpness
- Excellently clean, precise and force-free incisions
- Protection mechanism for storage of the blade inside the handle
- Color coded Titanium handles
- Four different cutting geometries: round, retrograde, wedge and lancet blade
- For further details see brochure no. C22402.



FD113D
Round blade,
gold-colored
7 facets
Length
205 mm, 8"




FD114D
Retro blade,
copper-colored
60°
Length
205 mm, 8"




FD115D
Wedge blade,
black-colored
45°
Length
205 mm, 8"



FD116D
Lancet blade,
bronze-colored
60°
Length
205 mm, 8"



SEM view of a diamond knife blade



SEM view of a common scalpel blade

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery NOIR[®] Brain Spatulas

Noir[®] Coated brain spatula

- Malleable metal spatula
- Conically tapered
- Smooth surface
- Atraumatic rounded edges
- Black Noir[®] surface coating to avoid light reflections
- Reusable
- Easy to reprocess



FF456B

S = 8 x 4 mm
Length:
200 mm, 8"



FF457B

M = 13 x 6 mm
Length:
200 mm, 8"



FF458B

L = 17 x 9 mm
Length:
200 mm, 8"



FF459B

XL = 21/11 mm
Length:
200 mm, 8"

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery
NOIR[®] coated KERRISON Detachable Bone Punches

Jaw position 130°, upbiting



Shaft length	Width	Footplate	NOIR [®] , detachable	Ejector	Jaw opening
180 mm, 7"	1.0 mm	thin	FK900B	-	8 mm
	1.5 mm	thin	FK911B	-	9 mm
	2.0 mm	thin	FK901B	✓	9 mm
	2.5 mm	thin	FK912B	✓	10 mm
	3.0 mm	thin	FK902B	✓	10 mm
200 mm, 8"	1.5 mm	standard	FK966B		9 mm
	2.0 mm	standard	FK913B		9 mm
	2.5 mm	standard	FK967B		10 mm
	3.0 mm	standard	FK914B	✓	10 mm



At a glance, large numbered jaw identification



Ejector - for the easy removal of punched-out material.




Numerical code - for reliable identification when assembling the two punch components.

MINOP[®] TEAM

Transcranial Endoscope Assisted Microneurosurgery
Bayonet-shaped KERRISON Punches

Jaw position 130°, upbiting



Shaft length	Width	Working length	Detachable	Jaw opening
240 mm, 9½"	2.0 mm	170 mm, 6¾	FF496R	10 mm
	3.0 mm	170 mm, 6¾	FF497R	10 mm
	4.0 mm	170 mm, 6¾	FF498R	10 mm
	5.0 mm	170 mm, 6¾	FF499R	10 mm




■ For more information about MINOP[®] TEAM please see our „Practical Atlas“ C29802.



■ For more information about KERRISON bone punches please see our brochure C84802.





TRANSNASAL NEUROENDOSCOPY

MINOP[®] TREND

TRansnasal ENDoscopic System





"When looking at recent publications on transsphenoidal surgery, it will be clear that **TR**anssphenoidal **END**oscopy is TREND-setting! However, this endoscopic technique is not in routine use everywhere and neurosurgeons are often reluctant to use it: One is often cautious about an endoscopic endonasal dissection because the permanent contamination of the endoscope with blood and nasal secretions hinders orientation. In addition, the para-endoscopic and biportal dissection is very unfamiliar requiring an unacceptably steep learning curve.

Nevertheless, endoscopic visualization and para-endoscopic dissection without using the surgical microscope offers several undisputable advantages. Advantages in visualization increases light intensity in the deep-seated surgical field and clearly displays patho-anatomical details. In addition, the extended viewing angle of endoscopes enables surgeons to observe hidden parts of the surgical field. The major benefit in surgical dissection is the unhindered approach to these clearly visible structures: Without using a nasal speculum, surgical manipulation is not impeded and the instruments are freely mobile. In addition, a pure endoscopic technique avoids the need

for rhinoseptal submucosal dissection providing a direct and quicker approach to the sphenoid sinus. This method avoids the need for postoperative nasal packing, thus causing less pain and discomfort after surgery, providing better nasal airflow and a shorter hospital stay.

Pre-conditions of transsphenoidal endoscopy are the basic endoscopic experience and anatomical studies in the laboratory; however, it is indispensable to use a dedicated endoscopic system to further shorten the learning phase. The endoscope for transsphenoidal skull base surgery must provide a brilliant image quality with true colors, high contrast and highly realistic images. This simplifies the differentiation between healthy or pathological structures. It is essential to have an effective cleaning function in order to free the endoscope lens from fog, blood or mucosal secretions. The endoscope must offer a highly ergonomic design and sufficient working length for extended approaches. For selected cases, it is also necessary to connect the endoscope to a navigation system or a holding device."

André Grotenhuis



André Grotenhuis
Nijmegen, Netherlands

MINOP® TREND

TRANSNASAL ENDOSCOPIC SYSTEM

MINOP® TREND Handle, Accessories and Storage

- Efficient suction, cleaning and irrigation function
- Highly ergonomic handle
- No irrigation pump needed
- Excellent image quality
- Different directions of view
- Big range of suitable instruments available



FH615

Ergonomic Handle with irrigation button
for FH610R and FH611R




RT099R

Adapter for Aesculap holding arm

FH605SU

Single use suction and irrigation tube, sterile packed,
Length 4.5 m, 2 puncture needles, for MINOP® TREND handle
FH615, Sales unit: PAK = Package of 10 tubes



FF357R

Storage tray with silicone padding and lid for all MINOP®
TREND components

only for reprocessing, not for transportation/shipment
(instruments not included)

Dimensions: (L/W/H) 410 x 257 x 64 mm



3/4 Container (basic version)

for storage racks FF357R, consisting of:

JK740

Container body 3/4 for FF357R without base perforation

Outside/Inside dimensions with lid:

L/W/H 470 x 285 x 112mm

L/W/H 421 x 258 x 75mm

JK789

Inner lid 3/4, blue

"The view through the operating microscope allows a purely coaxial visualisation in transsphenoidal surgery: laterally located structures are concealed behind the nasal speculum. Blind tumor removal involves a higher risk of iatrogenic damage to neurovascular structures and a possible increase in tumor remnants. With the use of the MINOP TREND endoscope for transnasal procedures, these laterally located parts of the field are directly visible and therefore surgically better approachable. In the past 15 years of endoscopic transnasal surgery, the use of endoscopes has proven to be not only indispensable but rather mandatory for a safe and effective transnasal surgery in the sellar and parasellar region."

André Grotenhuis, Nijmegen, Netherlands



MINOP[®] TREND

TRANSNASAL ENDOSCOPIC SYSTEM

MINOP[®] TREND Trocars and Endoscopes

FH610R

Suction and irrigation trocar

for 0° endoscope PE487A

Diameter: 4.5 / 6.0 mm

Working length: 120 mm



FH611R

Suction and irrigation trocar

for 30° endoscope PE507A

Diameter: 4.5 / 6.0 mm

Working length: 120 mm




PE487A

Endoscope

Direction of view: 0° (green ring)

Shaft diameter 4.0 mm

Autoclavable



PE507A

Endoscope

Direction of view: 30° (red ring)

Shaft diameter 4.0 mm

Autoclavable



"No other system that I have used combines as many helpful features in a single instrument". The lens cleaning is rapid and conveniently controlled with a button, instead of a pedal. The suction is effective. The ability to rotate the scope easily and quickly within the handle improves angled viewing. Overall, these features make the MINOP TREND an asset for endonasal surgery."

Jeremy Greenlee, Iowa City, USA

MINOP[®] TREND

TRansnasal ENDoscopic System - TREND – Currettes and Dissectors



FA041R-FA068R
Working length:
130 mm, 5 1/8"
Total length:
280 mm, 11"

diam. 6.5 mm		diam. 6.5 mm			
					
NICOLA FA041R	NICOLA FA042R	HARDY FA043R	HARDY FA044R		
Curette 45° vertical angled long neck	Curette 45° horizontal angled short neck	Enucleator left cutting	Enucleator right cutting		
diam. 4.0 mm		diam. 4.0 mm			
					
HARDY FA045R	HARDY FA046R	HARDY FA047R	HARDY FA060R		
Curette 90° left angled long neck	Curette 90° left angled short neck	Curette 90° right angled long neck	Curette 90° right angled short neck		

diam. 4.0 mm



HARDY
FA061R

Curette
45° horizontal,
left angled.
short neck

diam. 4.0 mm



HARDY
FA062R

Curette
45° horizontal,
right angled
short neck

diam. 6.0 mm



HARDY
FA063R

Curette
90° left angled
long neck

diam. 6.0 mm



HARDY
FA064R

Curette
90° left angled
short neck

diam. 6.0 mm



HARDY
FA065R

Curette
90° right angled
long neck

diam. 6.0 mm



HARDY
FA066R

Curette
90° right angled
short neck

diam. 1.7 mm



REULEN-
LANDOLT
FA067R

Micro Hook

diam. 2.0 mm



REULEN-
LANDOLT
FA068R

Dissector
blunt



MINOP[®] TREND

TRansnasal ENDoscopic System - TREND – Currettes and Dissectors

Straight design with ergonomic grasping part and semi-sharp tips



	diam. 6.5 mm	diam. 6.5 mm		diam. 4.0 mm	diam. 4.0 mm
					
	NICOLA FA030R	NICOLA FA031R	HARDY FA032R	HARDY FA033R	HARDY FA034R
	Curette 45° vertical angled, long neck	Curette 45° horizontal angled, short neck	Enucleator left cutting	Enucleator right cutting	Curette 90° angled long neck
					
					HARDY FA035R
					Curette 90° angled short neck
		diam. 4.0 mm	diam. 6.0 mm	diam. 6.0 mm	diam. 1.7 mm
					
		HARDY FA036R	HARDY FA037R	HARDY FA038R	LANDOLT-REULEN FA039R
		Curette 45° angled short neck	Curette 90° angled long neck	Curette 90° angled short neck	Micro Hook
					
					LANDOLT-REULEN FA040R
					Dissector blunt

FA030R-FA040R

Working length:
140 mm, 5 1/2"

Total length:
265 mm, 10 1/2"



MINOP[®] TREND

TRansnasal ENDoscopic System - Nasal Specula



1/2

COTTLE

OK105R-OK108R

OK090R

with aseptic joint, set-screw,
with extra thin blades
140 mm, 5 1/2"



OK105R 1/1




OK106R 1/1




OK107R 1/1



OK108R 1/1



OK090R 1/1



KILLIAN

- OK081R** Fig. 1
- OK082R** Fig. 2
- OK083R** Fig. 3
- OK084R** Fig. 4

with aseptic joint
145 mm, 5 3/4"

- OK091R** Fig. 1
- OK092R** Fig. 2
- OK093R** Fig. 3
- OK094R** Fig. 4

with screw joint
140 mm, 5 1/2"




Fig.1 1/1




Fig.2 1/1



Fig.3 1/1




Fig.4 1/1


MINOP[®] TREND

TRansnasal ENDoscopic System - Pituitary Instruments / Sinus Punches

FA076R

Antrum punch for removal of posterior nasal septum,

Rotating sheath 360°,
Working length: 120 mm, 4 3/4"



OK602R-OK609R

Sinus Punches




		 2/1 6 x 1.5 mm	 2/1 8 x 3 mm	 2/1 11.5 x 3.5 mm
	straight	OK608R forward through cutting	MACKAY-GRUNEWALD OK602R forward through cutting	MACKAY-GRUNEWALD OK603R forward through cutting
	45° upwards angled	OK609R forward through cutting	MACKAY-GRUNEWALD OK606R forward through cutting	MACKAY-GRUNEWALD OK607R forward through cutting

MINOP[®] TREND

TRansnasal ENDoscopic System - Antrum Grasping Forceps

OK680R

jaw opening backwards,
curved downwards



OK681R

jaw opening backwards,
curved upwards



OK682R

jaw opening backwards,
curved to right




OK683R

jaw opening backwards,
curved to left



OK684R

jaw opening backwards,
jaw 360° rotatable




MINOP[®] TREND

TRansnasal ENDoscopic System - Nasal Forceps

LANDOLT

FF345R


Tumor grasping forceps,
blunt, straight
Diam. 9.0 mm



TAKAHASHI

OK525R


Rongeur, straight



MINOP[®] TREND

TRansnasal ENDoscopic System - Nasal Forceps


Angled positions for rongeurs



BLAKESLEY-WILDE OK505R-OK509R

Ethmoidal forceps, straight


OK505R	3.0 mm		1/1
OK506R	3.6 mm		1/1
OK507R	4.2 mm		1/1
OK508R	4.8 mm		1/1
OK509R	5.6 mm		1/1



BLAKESLEY-WILDE OK520R-OK522R

Ethmoidal forceps, upwards curved, 140°

OK520R	3.6 mm		1/1
OK521R	4.2 mm		1/1
OK522R	4.8 mm		1/1



■ For more information about instruments for Functional Endoscopic Sinus Surgery (FESS), please ask your local Aesculap sales representative or see brochure no. C87511.

MINOP[®] TREND

TRansnasal ENDoscopic System - Nasal Scissors

OK560R - OK562R

Nasal scissors

OK560R


straight, blades serrated

OK561R

left curved, blades serrated

OK562R



right curved, blades serrated



CASPAR

FD228R

Micro scissors, curved
rotatable 360°,
detachable for cleaning



MINOP[®] TREND

TRansnasal ENDoscopic System - Pituitary Scissors

165 mm, 6 1/2"



FAHLBUSCH

FD220R

Micro scissors, extra delicate pattern, curved on flat, horizontal cutting



NICOLA

FD222R

Forceps, scoop-shaped, diam. 2.5 mm



YASARGIL-NICOLA

FD224R

Grasping forceps with long conical jaw



NICOLA

FD226R

Micro scissors, straight, diam. 2.5 mm




1/2

FD220R-FD226R

extra delicate tubular shaft scissors and grasping instruments for pituitary & skull base surgery

115 mm, 4 1/2"



CASPAR

FD228R

Micro scissors, curved rotatable sheath 360°

1/1



1/2





"Essential part of the endoscopic transnasal surgery is the nasal dissection, using special pituitary instruments. Goal is the maximum exploration of the target area, but also minimally invasive nasal traumatization, thus avoiding mucosal lacerations and unnecessary bony fractures. This influences patients postoperative quality of life enormously."

André Grotenhuis, Nijmegen, Netherlands

MINOP[®] TREND

TRansnasal ENDoscopic System - Pituitary Scissors





FA069R

^{1/1} straight



FA070R

^{1/1} right curved



FA071R

^{1/1} left curved

FA069R-FA071R
Micro Forceps


MINOP[®] TREND

TRansnasal ENDoscopic System - Curved Micro Suction Instruments

Suction cannulas

Curved suction instruments

FUKUSHIMA DESIGN




	Outer diameter	Inner diameter	Angled tip	Working length	Total length
GF431R	2.7 mm	2.0 mm	Right angled tip	135 mm, 5 1/4"	200 mm, 8"
GF432R	2.7 mm	2.0 mm	Left angled tip	135 mm, 5 1/4"	200 mm, 8"



MINOP[®] TREND

TRansnasal ENDoscopic System - Micro Suction Instruments



			○S	○○M	○○○L	○○○○XL
Working length			100 mm	115 mm	140 mm	165 mm
Total length			185 mm	200 mm	225 mm	250 mm
Outer diameter (3Fr = 1 mm)	3 Fr.	○	GF240R	GF250R	GF260R	GF270R
	4 Fr.	○	GF241R	GF251R	GF261R	GF271R
	5 Fr.	○	GF242R	GF252R	GF262R	GF272R
	6 Fr.	○	GF243R	GF253R	GF263R	GF273R
	7 Fr.	○	GF244R	GF254R	GF264R	GF274R
	8 Fr.	○	GF245R	GF255R	GF265R	GF275R
	9 Fr.	○	GF246R	GF256R	GF266R	GF276R
	10 Fr.	○	GF247R	GF257R	GF267R	GF277R
	12 Fr.	○	GF248R	GF258R	GF268R	GF278R

MINOP[®] TREND

TRansnasal ENDoscopic System - Nasal Forceps


GK801R

Bipolar coagulation forceps
with slender jaws and higher
spring tension

Total length 255 mm, 10"

Working length 135 mm, 5 1/4"

Aesculap tab connector — 




GK800R

T-coagulation forceps
with blunt, t-shaped tips

Total length 255 mm, 10"

Working length 135 mm, 5 1/4"

Aesculap tab connector — 



LANDOLT


GK560R

Coagulation forceps for
hypophysectomy, 90°

LANDOLT

GK580R

Coagulation forceps for
hypophysectomy, 120°



OF601R

Sickle knife, sharp tip
190 mm, 7 1/2"



BN175R

Frontal sinus ostium seeker,
double ended, curved
220 mm, 8 3/4"



FM158R


Micro tissue grasping forceps,
bayonet-shaped, straight tip
Total length 245 mm, 9 5/8"
Working length 120 mm, 4 3/4"



FM156R  Jaw 0.5 mm

FM157R  Jaw 0.9 mm

Micro Tissue grasping forceps
bayonet-shaped, straight tip
Working length 120 mm, 4 3/4"
Total length 240 mm, 9 1/2"



MINOP[®] TREND

TRansnasal ENDoscopic System - KERRISON Bone Punches

Jaw position 130°, upbiting



Shaft length	Width	Footplate	Detachable	Ejector	NOIR [®] , detachable	Ejector	Jaw opening
180 mm, 7"	1.0 mm	thin	FK906R	-	FK906B	-	8 mm
	1.5 mm	thin	FK923R	-	FK923B	-	9 mm
	2.0 mm	thin	FK907R	✓	FK907B	✓	9 mm
	2.5 mm	thin	FK924R	✓	FK924B	✓	10 mm
	3.0 mm	thin	FK908R	✓	FK908B	✓	10 mm
	4.0 mm	thin	FK909R	✓	FK909B	✓	12 mm

Jaw position 130°, downbiting




Shaft length	Width	Footplate	Detachable	Ejector	Jaw opening
180 mm, 7"	1.0 mm	thin	FK936R	-	8 mm
	1.5 mm	thin	FK937R	-	9 mm
	2.0 mm	thin	FK938R	✓	9 mm



MINOP[®] TREND

TRansnasal ENDoscopic System - KERRISON Bayonet Bone Punches

Jaw position 130°, upbiting



Shaft length	Width	Working length	Detachable	Jaw opening
240 mm, 9½"	2.0 mm	170 mm, 6¾"	FF496R	10 mm
	3.0 mm	170 mm, 6¾"	FF497R	10 mm
	4.0 mm	170 mm, 6¾"	FF498R	10 mm
	5.0 mm	170 mm, 6¾"	FF499R	10 mm



■ For more information about MINOP[®] TREND please see our „Practical Atlas“ C26402.

AESCULAP NEUROSURGERY

M-TRAC® – Mechanical Holding Device

M-TRAC® – Flexible holding device with mechanical fixation

- Total length: 107 cm
- Length of fixation bar: 46 cm
- Diameter of fixation bar: 20 mm
- Total weight: 0.7 kg
- Holding force: 4 kg
- Easy mechanical fixation by clamping handle
- Small, flexible joints for fine positioning
- Autoclavable 134°C, 5 minutes
- Full range of accessories/adapters for connecting Aesculap endoscopes, trocars and instruments
- Holding arm fits into regular Standard 1/1 Container, see brochure no. C40402



FF280R

Flexible fixation element with ball joint suitable for RT040R and FF168R



RT090R

Flexible fixation element with sprocket suitable for RT040R and FF168R



FF151R

Rigid fixation element suitable for RT040R and FF168R



- For further details see brochure no. C26911



AESCULAP NEUROSURGERY

UNITRAC® – Pneumatic Holding Device

UNITRAC® – Pneumatic holding arm

- Single handed use
- Full range of accessories / adapters for connecting Aesculap endoscopes, trocars and instruments
- Integrated safety systems prevent collapse of holding arm if OR compressed air supply is interrupted
- Direct connection to OR compressed air supply
- Diameter of fixation bar: 20 mm
- To be used with JG901



JG901

Sterile drape for coverage of the Unitrac® arms, single-use product,

Sale unit:
PAK = Package of 50 pieces



RT020R

Quick connect adapter for use with sterile drape JG901 allows the change of instruments after draping with JG901



RT043R

CO₂ cartridge adapter for use of UNITRAC®, independent from compressed air sources



RT044SU

CO₂ air cartridge, single use

Sale unit:
PAK = Package of 10 pieces



- For further details see brochure no. C47411

AESCLAP NEUROSURGERY

Adapters for UNITRAC® and M-TRAC®

RT046P

Universal holder

for endoscopes and trocars with diam. 3.0-7.5 mm, consisting of: RT081R and RT055P



RT081R

Adapter

for universal insert RT055P



RT055P

Universal insert (Spare Part)

for endoscopes and trocars with diam. 3.0-7.5 mm



RT099R

Adapter

for fixation of MINOP® TREND handle, FH615



RT079R

Adapter

for fixation of angled endoscopes PE486A, PE506A, PE526A



RT068R

Adapter

for fixation of MINOP® InVent trocar, FH620R



RT079205

Silicone bit

for RT079R



	MINOP® FF397R FF398R FF399R	Paediscopes PA010A	MINOP® InVent FH620R	Angled scopes PE486A PE506A PE526A	MINOP® TREND FH615	MINOP® TR FH601R
RT046P	●	●		●		●
RT099R					●	
RT079R				●		
RT068R			●			

AESCULAP NEUROSURGERY

Holding Devices - NEUROPILOT® – Fine-positioning for UNITRAC® and M-TRAC®

NEUROPILOT® for IntraVentricular and endoscope-assisted indications. NEUROPILOT® is a new, unique steering device for neuroendoscopes. After positioning the neuroendoscope in situ, finest corrections or adjustments are necessary, to receive the optimal endoscopic image. With traditional holding devices, only rough positioning is possible; a precise and fine steering of the neuro-endoscope can be compromised.

NEUROPILOT® offers a number of unique advantages:

- Proper fixation of the neuro-endoscope in the NEUROPILOT® and the holding device
- Precise steering of the neuro-endoscope by three screws in the three-dimensional space
- Accurate manoeuvring of the neuro-endoscope by defined movements in the sub-millimeter area



"In pure intraventricular neuroendoscopy, a micro-steering device can be extremely useful. If the precision and adjustment of a holding arm is not enough, the Neuropilot closes this gap. Additionally, in cases where both hands are needed for instrumentation the Neuropilot is of great help.

The Aesculap Neuropilot is the only system on the market providing finest correction of your endoscope in a three-dimensional space inside the ventricular compartments."

Peter Nakaji, Phoenix, USA



RT060R

NEUROPILOT®

for intraventricular and
endoscope-assisted indications



RT061R

Insert for angled endoscopes
PE486A - PE526A with diam. 4.0 mm



RT065R

Insert for MINOP® trocar FF399R
with diam. 6.0 mm



RT063R

Insert for MINOP® trocar FF397R
with diam. 3.2 mm



RT066R

Insert for PaediScope® PF010A
with diam. 3.0 mm



RT064R

Insert for MINOP® trocar FF398R
with diam. 4.6 mm



AESCULAP NEUROSURGERY

Holding Devices - NEUROPILOT® – Fine-positioning for UNITRAC® and M-TRAC®

	MINOP® FF397R	MINOP® FF398R	MINOP® FF399R	Paediscope PA010A	Angled scopes PE486A PE506A PE526A	MINOP® TR FH601R
RT060R	●	●	●	●	●	●
RT061R					●	
RT063R	●					
RT064R		●				●
RT065R			●			
RT066R				●		

AESCULAP NEUROSURGERY

Operating Manuals



Aesculap® MINOP®

Intraventricular Neuroendoscopy: A Practical Atlas

Mark M. Souweidane, M.D., F.A.C.S, F.A.A.P.

C29202

Helsinki Microneurosurgery App

App Store > Search > B. Braun AG

B|BRAUN SHARING EXPERTISE



iTunes



Google Play

Aesculap Neuroendoscopy App

App Store > Search > B. Braun AG

B|BRAUN SHARING EXPERTISE



iTunes



Google Play

AESCLAP NEUROSURGERY

Visual Equipment Examples for Neuroendoscopy

Full HD Camera, Xenon Light Source, Documentation System and Touch Screen

PV463

3-Chip Full HD camera head
with fixed focus coupler

PV462

3-Chip Full HD camera head
with zoom-coupler

PV884

Camera holder

PV460

Full HD Cameracontrol unit

PV880

"Metro Junior" Endoscopy cart
835 x 1580 x 750 mm (w x h x d)

PV881

"Metro Junior" Endoscopy cart
with integrated Isolation Transformer
835 x 1580 x 750 mm (w x h x d)



PV956

24" Full HD touch panel
display

PV909

Monitor stand
for PV959 and PV956

OP940

LED light source

PV840

EDDY DVD
Digital Documentation System

JG904

Sterile Camera drape,
disposable, ring design,
package of 25



AESCLAP NEUROSURGERY

Visual Equipment Examples for Neuroendoscopy

Full HD CMOS Camera, LED Light Source and Flat Screen

PV473

Full HD CMOS camera head with fixed focus coupler

PV472

Full HD CMOS camera head with zoom-coupler

PV884

Camera holder

PV460

Full HD Cameracontrol unit

PV959

26" Full HD touch panel display

PV909

Monitor stand for PV959 and PV956

OP940

LED light source

PV880

"Metro Junior" Endoscopy cart
835 x 1580 x 750 mm (w x h x d)

PV881

"Metro Junior" Endoscopy cart with integrated Isolation Transformer
835 x 1580 x 750 mm (w x h x d)



OP923

Full HD Light cable, autoclavable, diam. 4.8 mm, length 250 cm



■ For further details see brochure no. C46702





POWER SYSTEMS FOR NEUROSURGERY

AESCULAP NEUROSURGERY

Power Systems - ELAN® 4 electro – Electric Highspeed Power System

GA800

ELAN® 4 electro control unit



GA806

ELAN® 4 electro motor cable



GA808

ELAN® 4 electro foot control



GA861

**ELAN® 4 electro 1-ring handpiece
L4**



GA862

**ELAN® 4 electro 1-ring handpiece
L7**



GA863

**ELAN® 4 electro 1-ring handpiece
L10**



GA864

**ELAN® 4 electro 1-ring handpiece
L13**



AESCALAP NEUROSURGERY

Power Systems - ELAN® 4 air - Pneumatic Highspeed Power System

GA702R

ELAN® 4 air wall hose
3 m



GA710R

WALL ADAPTER
Aesculap-Draeger



GA712R

WALL ADAPTER Schrader



GA713R

WALL ADAPTER DIN



GA715R

WALL ADAPTER Synthes

GA708

ELAN® 4 air foot control



GA705R

ELAN® 4 air motor hose
3 m for foot control



GA706R

ELAN® 4 air motor hose
3 m with hand control



GA707R

ELAN® 4 air motor hose
5 m with hand control

GA761

ELAN® 4 air 1-ring handpiece L4



GA762

ELAN® 4 air 1-ring handpiece L7



GA763

ELAN® 4 air 1-ring handpiece L10




GA764


ELAN® 4 air 1-ring handpiece L13



AESCLAP NEUROSURGERY


Power Systems - ELAN® 4 tools for 1-ring handpieces

Ø	0.6 mm	0.8 mm	1.0 mm	1.0 mm	1.4 mm	1.4 mm	1.8 mm	1.8 mm	2.3 mm	2.3 mm	2.7 mm	
												 sterile
I	GP111R	GP112R	GP113R	GP114R	GP115R	GP116R	GP117R	GP118R	GP120R	GP121R	GP122R	Rosen


Ø	2.7 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm					
												 sterile
I	GP123R	GP124R	GP125R	GP126R	GP127R	GP128R	GP129R					Rosen

Ø	3.0 mm	4.0 mm	5.0 mm	6.0 mm	
Soft cut 					 sterile
I	GP133R	GP134R	GP135R	GP136R	Rosen

Ø	4.0 mm	5.0 mm	6.0 mm	7.0 mm	
					 sterile
I	GP184R	GP185R	GP186R	GP187R	Twin-Cut burr



Ø	0.6 mm	0.8 mm	1.0 mm	1.0 mm	1.4 mm	1.4 mm	1.8 mm	1.8 mm	2.3 mm	2.3 mm	2.7 mm	
												 sterile
I	GP141R	GP142R	GP143R	GP144R	GP145R	GP146R	GP147R	GP148R	GP149R	GP150R	GP151R	Diamond burr

Ø	2.7 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm					
												 sterile
I	GP152R	GP153R	GP154R	GP155R	GP156R	GP157R	GP158R					Diamond burr


AESCLAP NEUROSURGERY

Power Systems - ELAN[®] 4 tools for 1-ring handpieces

Ø	2.0 mm	2.0 mm	2.3 mm	2.3 mm	3.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm			
												
I	GP161R	GP162R	GP163R	GP164R	GP165R	GP166R	GP168R	GP169R	GP170R			
												 sterile Diamond burr coarse

Ø	3.0 mm	4.0 mm	5.0 mm	6.0 mm	7.0 mm							
												
I	GP173R	GP174R	GP175R	GP176R	GP177R							
												 sterile Diamond burr extra coarse





Ø	1.5 mm	2.0 mm	2.5 mm	3.0 mm	3.0 mm
US Type					
I	GP201R	GP202R	GP203R	GP204R	GP205R

sterile
Neuro cutter

Ø	2.0 mm	2.5 mm	3.0 mm
Aesculap Type			
I	GP208R	GP209R	GP210R

sterile
Neuro cutter

Ø	1.5 mm	1.5 mm	2.0 mm	2.0 mm	2.5 mm	2.5 mm	3.0 mm	3.0 mm	4.0 mm	4.0 mm		
I	GP211R	GP212R	GP213R	GP214R	GP215R	GP216R	GP217R	GP218R	GP219R	GP220R		

sterile
Neuro cutter diamond

Ø	4.0 mm	5.0 mm	6.0 mm
I	GP224R	GP225R	GP226R

sterile
Barrel burr

Ø	4.0 mm	5.0 mm	6.0 mm
Soft cut			
I	GP228R	GP229R	GP230R


sterile
Barrel burr

Ø	4.0 mm	5.0 mm	6.0 mm
I	GP234R	GP235R	GP236R

sterile
Cone burr

AESCULAP NEUROSURGERY

Power Systems - ELAN® 4 tools for 1-ring handpieces



Ø	3.0 mm	4.0 mm	5.0 mm	6.0 mm
Midas Type				
	GP193R	GP194R	GP195R	GP196R

sterile

Acorn burr

Ø	4.0 mm	5.5 mm
	GP232R	GP233R

sterile

Oval burr

Ø	4.0 mm
	GP238R

sterile

Reverse taper burr coarse diamond

Ø	1.0 mm	2.0 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm
Tungsten Carbide						
	GP113TC	GP119TC	GP124TC	GP126TC	GP128TC	GP129TC

sterile

Rosen

Ø	3.0 mm
Tungsten Carbide	
	GP210TC

sterile

Neuro cutter

Ø	2.3 mm	1.8 mm	2.1 mm	2.3 mm
	GP240R	GP241R	GP242R	GP243R

sterile

Lindemann

Ø	1.5 mm	2.0 mm
	GP188R	GP189R

sterile

Twist drill

Ø	1.0 mm
	GP190R

sterile

Pin cutter



Power Systems - ELAN® 4 craniotome- / multifunction handpieces, attachments and tools

GA849

ELAN® 4 electro craniotome and multifunction handpiece 2-ring



GA749

ELAN® 4 electro craniotome and multifunction handpiece 2-ring



GB945R

Holding sleeve



GB941R

Fixed dura guard
PAEDIATRIC



GB942R

Fixed dura guard
STANDARD



GB943R

Fixed dura guard
LONG




GB947R

Turnable dura guard
STANDARD

AESCULAP NEUROSURGERY

Power Systems - ELAN® 4 tools for craniotome and multifunction handpieces




∅	2.3 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm	7.0 mm	8.0 mm			
II	GP301R	GP302R	GP303R	GP304R	GP305R	GP306R	GP307R	GP308R	GP309R			1 sterile Rosen

∅	2.3 mm	3.0 mm	3.5 mm	4.0 mm	4.5 mm	5.0 mm	6.0 mm					
II	GP311R	GP312R	GP313R	GP314R	GP315R	GP316R	GP317R					1 sterile Diamond burr


∅	2.3 mm	3.0 mm	4.0 mm	5.0 mm	6.0 mm							
II	GP321R	GP322R	GP323R	GP324R	GP325R							1 sterile Diamond burr coarse

∅	4.0 mm	5.0 mm	6.0 mm									
II	GP328R	GP329R	GP330R									1 sterile Diamond burr extra coarse




 sterile

Twist drill




 sterile

Twist drill




 sterile

Acorn burr




 sterile

Twin-Cut burr



 sterile

Craniotome cutter - spiral type



 sterile

Craniotome cutter - straight type

- The whole burr range for 2-ring handpieces can be found in the Burrs & Blades catalogue 017599.

AESCULAP NEUROSURGERY

Power Systems - ELAN® 4 perforator drivers and tools & ELAN® 4 lowspeed motors

GA822

ELAN® 4 electro perforator driver



GA722

ELAN® 4 air perforator driver



Skull perforators



Ø/mm


GB300R	6/9
GB302R	9/12
GB304R	12/15



Ø/mm

TE561	6/9
TE562	9/12
TE563	12/15

Hudson



GA824

ELAN® 4 electro lowspeed motor



GA724

ELAN® 4 air lowspeed motor



GA836

ELAN® 4 electro micro sagittal saw




GA736

ELAN® 4 air micro sagittal saw




GP491R	14 mm	13 mm	0.3 mm	0.3 mm
GP492R	15 mm	5 mm	0.3 mm	0.3 mm
GP493R	20 mm	5 mm	0.3 mm	0.3 mm
GP494R	20 mm	10 mm	0.3 mm	0.3 mm
GP495R	20 mm	15 mm	0.3 mm	0.3 mm
GP496R	25 mm	5 mm	0.3 mm	0.3 mm
GP497R	25 mm	12 mm	0.3 mm	0.3 mm

GP491R




GP495R



GP492R




GP496R




GP493R



GP497R



GP494R



AESCU LAP NEUROSURGERY

Power Systems - ELAN® 4 micro reciprocating saws and saw blades

GA837

ELAN® 4 electro micro reciprocating saw




GA737

ELAN® 4 air micro reciprocating saw




GP542R	13 mm	0.3 mm	0.3 mm
GP543R	20 mm	0.3 mm	0.3 mm
GP544R	13 mm	0.3 mm	0.3 mm
GP545R	20 mm	0.3 mm	0.3 mm


GP542R




GP543R



GP544R



GP545R



Power Systems - ELAN 4 Accessories

TE780

POWER CORD Europe, 1.5 m

TE730

POWER CORD Europe, 5.0 m

TE734

POWER CORD UK, 5.0 m

TE734

POWER CORD USA, Canada, Japan,
3.5 m



TE734

ELAN® 4 electro disposable tube set



GA259SU

ELAN® 4 spray nozzle
for craniotome attachment

GA261SU

ELAN® 4 spray nozzle
for 1-ring handpiece L4

GA262SU

ELAN® 4 spray nozzle
for 1-ring handpiece L7

GA263SU

ELAN® 4 spray nozzle
for 1-ring handpiece L10

GA264SU

ELAN® 4 spray nozzle
for 1-ring handpiece L13



GA258SU

ELAN® 4 spray nozzle
for saws




AESCLAP NEUROSURGERY

Power Systems - ELAN 4 cleaning and maintenance



GB691R

ELAN® 4 air rinsing device
for mechanical cleaning
not suitable for autoclaving



GB692R

ELAN® 4 electro rinsing device
for mechanical cleaning
not suitable for autoclaving



GB697R

ELAN® 4 air rinsing adaptor
for manual cleaning



GB698R

ELAN® 4 electro rinsing adaptor
for manual cleaning



GB600860

ELAN® 4 electro spray adaptor



GB600

STERILIT® Power Systems oil spray



GB600850

ELAN® 4 air spray adaptor

Main Catalogues

084002
Power Systems



017599
Burrs & Blades



022702
Power Systems – Old Systems*



ELAN 4

071602
ELAN® 4



071802
ELAN® 4 Burr Compendium*



Elan 4 Reprocessing

071511
Rinsing Device



012002
Manual Cleaning




038902
Acculan 3Ti Electro



AESCULAP ACADEMY

Forum for contemporary medicine.



Technical developments, new treatment methods, hospital management – the requirements placed on today's medical professionals are varied, which is why high-quality continuing professional development is more important than ever. This is precisely what Aesculap Academy stands for.

Its aim is clear: The Aesculap Academy wants to keep medical and specialist staff in hospitals and practices fit for the future. Founded in 1995 under the B. Braun Group, the Aesculap Academy is seen today around the world as an important forum for medical training and further professional education. It works with an interdisciplinary, independent and international approach and it strives for long-term partnerships. Thanks to tailored and interrelated modules, the participants can continually develop and build up their knowledge and skills throughout their career. This means they are always as well prepared as possible for their daily work and future tasks.

Take part in one of our international Neuroendoscopy courses.

For detailed information and registration please visit our website „www.aesculap-neuro.com“ or „www.aesculap-academy.com“ or contact your local B. Braun Aesculap representative.

AESCULAP ACADEMY

Am Aesculap-Platz 78532 Tuttlingen Phone +49 7461 95-2001
www.aesculap-academy.com



Tuttlingen | Berlin | Bochum

AESCLAP NEUROSURGERY

Literature

M. M. Souweidane, P. F. Morgenstern, S. Kang et al.

Endoscopic Third Ventriculostomy in Patients with a Diminished Prepontine Interval

Journal of Neurosurgery: Pediatrics, Vol. 5, 250-254, March 2010

O. Sacko, S. Boetto, V. Lauwers-Cances, et al.

Endoscopic Third Ventriculostomy: Outcome Analysis in 368 Procedures

Journal of Neurosurgery: Pediatrics, Vol. 5, 68-74, January 2010

N. Luther, W. R. Stetler Jr., Ira. J. Dunkel, et al.

Subarachnoid Dissemination of Intraventricular Tumors Following Simultaneous Endoscopic Biopsy and Third Ventriculostomy

Journal of Neurosurgery: Pediatrics, Vol. 5, 61-67, January 2010

G. P. Lekovic, J. F. Kerrigan, S. Wait, et al.

In Situ Single-Unit Recording of Hypothalamic Hamartomas Under Endoscopic Direct Visualization

Neurosurgery, Vol. 65, Nr. 6, E1195-E1196, December 2009

B. D. Kollroy, F. A. Ponce, Scott D. Wait, et al.

Endoscopic Intraventricular Biopsy of Infundibular Langerhans Cell Histiocytosis: Case Report

Neurosurgery, Vol. 65, Nr. 1, E214-E215, July 2009

P. Pillai, M. N. Baig, Ch. S. Karas, et al.

Endoscopic Image-Guided Transoral Approach to the Craniovertebral Junction: An Anatomic Study Comparing Surgical Exposure and Surgical Freedom Obtained with the Endoscope and the Operating Microscope

Neurosurgery, Operative Neurosurgery 2, Vol. 64, ONS437-ONS444, May 2009

N. J. Hopf, A. Stadie, R. Reisch, et al.

Surgical Management of Bilateral Middle Cerebral Artery Aneurysms via a Unilateral Supraorbital Key-Hole Craniotomy

Minimally Invasive Neurosurgery, Vol. 52, 126-131, 2009

R. Reisch, A. Stadie, R. Kockro, et al.

The Minimally Invasive Supraorbital Subfrontal Key-Hole Approach for Surgical Treatment of Temporomesial Lesions of the Dominant Hemisphere

Minimally Invasive Neurosurgery, Vol. 52, 163-169, 2009

J. Leonardo, R. A. Hanel, W. Grand

Endoscopic Tracking of a Ventricular Catheter for Entry into the Lateral Ventricle: Technical Note

Minimally Invasive Neurosurgery, Vol. 52, 287-289, 2009

A. T. Stadie, R. Reisch, R. A. Kockro, et al.

Minimally Invasive Cerebral Cavernoma Surgery Using Keyhole Approaches – Solutions for Technique-Related Limitations

Minimally Invasive Neurosurgery, Vol. 52, 9-16, 2009

P. Pillai, M. Lubow, A. Ortega, et al.

Endoscopic Transconjunctival Surgical Approach to the Optic Nerve and Medial Intraconal Space: A Cadaver Study

Neurosurgery, Operative Neurosurgery 2, Vol. 63, OBS204-ONS209, October 2008

S. C. Froelich, K. M. Abdel Aziz, P. D. Cohen, et al.

Microsurgical and Endoscopic Anatomy of Lilliequist's Membrane: A Complex and Variable Structure of the Basal Cisterns

Neurosurgery, Operative Neurosurgery 1, Vol. 63, ONS1-ONS9, July 2008

B. C. Ong, P. A. Gore, M. B. Donnellan, et al.

Endoscopic Sublabial Transmaxillary Approach to the Rostral Middle Fossa

Neurosurgery, Operative Neurosurgery 1, Vol. 62, 30-37, March 2008

J. D. W. Greenlee, C. Teo, A. Ghahreman, et al.

Purely Endoscopic Resection of Colloid Cysts

Neurosurgery, Operative Neurosurgery 1, Vol. 62, ONS51-ONS56, March 2008

AESCLAP NEUROSURGERY

Literature

P. A. Gore, L. F. Gonzalez, H. L. Rekate, et al.

Endoscopic Supracerebellar Infratentorial Approach for Pineal Cyst Resection: Technical Case Report

Neurosurgery, Operative Neurosurgery 1, Vol. 62, March 2008

P. Cappabianca, G. Cinalli, M. Gangemi, et al.

Application of Neuroendoscopy to Intraventricular Lesions

Neurosurgery, Supplement, Vol. 62, No. 2, SHC575-SHC598, February 2008

J. P. Greenfield, L. Z. Leng, U. Chaudhry, et al.

Combined Simultaneous Endoscopic Transsphenoidal and Endoscopic Transventricular Resection of a Giant Pituitary Macroadenoma

Minimally Invasive Neurosurgery, Vol. 51, 306-309, 2008

P. Y. Hwang, C. Long Ho

Neuronavigation Using an Image-Guided Endoscopic Transnasal-Sphenoethmoidal Approach to Clival Chordomas

Neurosurgery, Operative Neurosurgery 2, Vol. 61, ONS212-ONS218, November 2007

B. Depreitere, N. Dasi, J. Rutka, et al.

Endoscopic Biopsy for Intraventricular Tumors in Children

Journal of Neurosurgery: Pediatrics, Vol. 106, 340-346, May 2007

M. Gangemi, F. Maiuri, G. Colella, et al.

Is Endoscopic Third Ventriculostomy an Internal Shunt Alone?

Minimally Invasive Neurosurgery, Vol. 50, 47-50, 2007

M. Husain, M. Rastogi, B. K. Ojha, et al.

Endoscopic Transoral Surgery for Craniovertebral Junction Anomalies

Journal of Neurosurgery: Spine, Vol. 5, 367-373, October 2006

C. Teo, P. Nakaji, R. J. Mobbs

Endoscope-Assisted Microvascular Decompression for Trigeminal Neuralgia: Technical Case Report

Neurosurgery, Operative Neurosurgery 4, Vol. 59, October 2006

A. Weyerbrock, T. Mainprize, J. T. Rutka

Endoscopic Fenestration of a Symptomatic Cavum Septum Pellucidum: Technical Case Report

Neurosurgery, Operative Neurosurgery 4, Vol. 59, October 2006

R. Moftakhar, M. S. Salamat, S. Sahin, et al.

Endoscopically-Assisted Resection of a Choroid Plexus Vascular Malformation Traversing the Cerebral Aqueduct: Technical Case Report

Neurosurgery, Operative Neurosurgery 1, Vol. 59, July 2006

J. van Beijnum, P. W. Hanlo, K. Sen Han, et al.

Navigated Laser-Assisted Endoscopic Fenestration of a Suprasellar Arachnoid Cyst in a 2-Year-Old Child with Bobble-Head Doll Syndrome

Journal of Neurosurgery: Pediatrics, Vol. 104, 348-351, May 2006

F. T. Mangano, D. D. Limbrick, J. R. Leonard, et al.

Simultaneous Image-Guided and Endoscopic Navigation without Rigid Cranial Fixation: Application in Infants: Technical Case Report

Neurosurgery, Operative Neurosurgery 2, Vol. 58, ONS-377-ONS-378, April 2006

G. P. Lekovic, L. F. Gonzalez, I. Feiz-Erfan, et al.

Endoscopic Resection of Hypothalamic Hamartoma using a Novel Variable Aspiration Tissue Resector

Neurosurgery, Operative Neurosurgery 1, vol. 58, ONS166-ONS169, February 2006

A. Morita, M. Shin, L. N. Sekhar, et al.

Endoscopic Microneurosurgery: Usefulness and Cost-Effectiveness in the Consecutive Experience of 210 Patients

Neurosurgery, Vol. 58, No. 2, 315-321, February 2006

A. A. Figaji, A. G. Fieggen, P. L. Semple, et al.

Intracranial Endoscopy

Samj Forum, Vol. 96, No. 1, 32-37, January 2006

A. Bussarsky, M. Marinov, V. Bussarsky, et al.

Virtual Simulation of Neuroendoscopic Procedures: Early Clinical Experience with Ventricular Lesions

Central European Neurosurgery, Vol. 67, 129-136, 2006

J. Zhao, Y. Wang, Y. Zhao, et al.

Neuroendoscope-Assisted Minimally Invasive Microsurgery for Clipping Intracranial Aneurysms

Minimally Invasive Neurosurgery, Vol. 49, 335-341, 2006

Y. Arakawa, K. Nakazawa, H. Kataoka, et al.

Microfiberscope Coaxial Technique in Neuroendoscopic Surgery

Minimally Invasive Neurosurgery, Vol. 49, 380-383, 2006

I. Gawish, R. Reisch, A. Perneczky

Endoscopic Aqueductoplasty through a Tailored Craniocervical Approach

Journal of Neurosurgery, Vol. 103, 778-782, November 2005

R. Reisch, A. Perneczky

Ten-Year Experience with the Supraorbital Subfrontal Approach through an Eyebrow Skin Incision

Neurosurgery, Operative Neurosurgery 4, Vol. 57, ONS242-ONS255, October 2005

M. M. Souweidane

Endoscopic Surgery for Intraventricular Brain Tumors in Patients without Hydrocephalus

Neurosurgery, Operative Neurosurgery 4, Vol. 57, ONS312-ONS318, October 2005

K. Schmidt, C. Coimbra

Endoscopic Treatment of Thalamic Neuroepithelial Cysts

Journal of Neurosurgery, Vol. 103, 342-346, August 2005

N. Luther, A. Cohen, M. M. Souweidane

Hemorrhagic Sequelae from Intracranial Neuroendoscopic Procedures for Intraventricular Tumors

Neurosurgical Focus, Vol. 19 (1), E9, 1-4, July 2005

P. D. Purdy, T. Fujimoto, R. E. Replogle, et al.

Percutaneous Intraspinal Navigation for Access to the Subarachnoid Space: Use of Another Natural Conduit for Neurosurgical Procedures

Neurosurgical Focus, Vol. 19 (1), E11, 1-5, July 2005

M. M. Souweidane

Endoscopic Management of Pediatric Brain Tumors

Neurosurgical Focus, Vol. 18 (6a), E1, June 2005

N. Luther, M. M. Souweidane

Neuroendoscopic Resection of Posterior Third Ventricular Ependymoma

Neurosurgical Focus, Vol. 18 (6a), E3, 1-2, June 2005

E. Nathal, J. L. Gomez-Amador

Anatomic and Surgical Basis of the Sphenoid Ridge Keyhole Approach for Cerebral Aneurysms

Neurosurgery, Operative Neurosurgery 1, Vol. 56, ONS178-ONS185, January 2005

J. Martin, C. Neal, I. Moores, et al.

Use of a Nitrogen Arm-Stabilized Endoscopic Microdriver in Neuroendoscopic Surgery

Minimally Invasive Neurosurgery, Vol. 48, 63-65, 2005

M.J. Fritsch, L. Dörner, S. Kienke, et al.

Hydrocephalus in children with posterior fossa tumors: The role of Endoscopic Third Ventriculostomy (ETV)

Journal of Neurosurgery: Pediatrics, Vol. 103: 40-42, 2005

AESCLAP NEUROSURGERY

Literature

M.J. Fritsch, S. Kienke, T. Ankermann, et al.

Endoscopic Third Ventriculostomy (ETV) in infants

Journal of Neurosurgery: Pediatrics, Vol. 103: 50-53, 2005

M. Taniguchi, A. Kato, T. Taki, et al.

Endoscope Assisted Removal of Jugular Foramen Schwannoma; Report of 3 Cases

Minimally Invasive Neurosurgery, Vol. 48, 365-368, 2005

H. Kinouchi, T. Yanagisawa, A. Suzuki, et al.

Simultaneous Microscopic and Endoscopic Monitoring During Surgery for Internal Carotid Artery Aneurysms

Journal of Neurosurgery, Vol. 101, 989-995, December 2004

J. C. Wang, L. Heier, M. M. Souweidane

Advances in the Endoscopic Management of Suprasellar Arachnoid Cysts in Children

Journal of Neurosurgery: Pediatrics, Vol. 100, 418-426, May 2004

M.J. Fritsch, S. Kienke, H.M. Mehdorn

Endoscopic aqueductoplasty: stent or not to stent?

Childs Nerv System, Vol. 20, 137-142, 2004

M.J. Fritsch, K.H. Manwaring, S. Kienke, et al.

Endoscopic treatment of isolated 4th ventricle in children

Neurosurgery, Vol. 55, 372-379, 2004

C. Trantakis, J. Helm, M. Keller, et al.

Third Ventriculostomy in Communicating Hydrocephalus in Adult Patients – The Role of Lumbar and Cranial Cerebrospinal Fluid Outflow Measurement

Minimally Invasive Neurosurgery, Vol. 48, 140-144, 2004

S. Wolfsberger, M.-T. Forster, M. Donat, et al.

Virtual Endoscopy is a Useful Device for Training and Preoperative Planning of Transsphenoidal Endoscopic Pituitary Surgery

Minimally Invasive Neurosurgery, Vol. 47, 214-220, 2004

Y. Lin, Y. Qiu

Microanatomy of Endoscope-Assisted Glabellar Nasal Keyhole Approach

Minimally Invasive Neurosurgery, Vol. 46, 155-160, 2003

T. G. Psarros, J. Krumerman, C. Coimbra

Endoscopic Management of Supratentorial Ventricular Neurocysticercosis: Case Series and Review of the Literature

Minimally Invasive Neurosurgery, Vol. 46, 331-334, 2003

M. A. Barajas, G. Ramirez-Guzmán, C. Rodríguez-Vazquez, et al.

Multimodal Management of Craniopharyngiomas: Neuroendoscopy, Microsurgery, and Radiosurgery

Journal of Neurosurgery (Supplement 5), Vol. 97, 607-609, December 2002

Z. Horváth, F. Vető, I. Balás, et al.

Biportal Endoscopic Removal of a Primary Intraventricular Hematoma: Case Report

Minimally Invasive Neurosurgery, Vol. 43, 4-8, 2000

A. Rieger, N. G. Rainov, M. Brucke, et al.

Endoscopic Third Ventriculostomy is the Treatment of Choice for Obstructive Hydrocephalus due to Pediatric Tumors

Minimally Invasive Neurosurgery, Vol. 43, 83-86, 2000

V. Rohde, J. M. Gilsbach

Anomalies and Variants of the Endoscopic Anatomy for Third Ventriculostomy

Minimally Invasive Neurosurgery, Vol. 43, 111-117, 2000

J. Paladino K. Rotim, D. Štimac, et al.

Endoscopic Third Ventriculostomy with Ultrasonic Contact Microprobe

Minimally Invasive Neurosurgery, Vol. 43, 132-134, 2000

T. Menovsky, J. A. Grotenhuis, J. de Vries, et al.

Endoscope-Assisted Supraorbital Craniotomy for Lesions of the Interpeduncular Fossa Technique and Application

Neurosurgery, Vol. 44, No. 1, 106-112, January 1999

P. Wieneke, T. Lutze

**Technologies for Microendoscopes of the Future:
The MINOP Project**

Minimally Invasive Therapy & Allied Technology, Vol. 7/3, 233-239, 1998

E. van Lindert, N. Hopf, A. Perneczky

**Endoscopic Treatment of Mesencephalic Ependymal Cysts:
Technical Case Report**

Neurosurgery, Vol. 43, No. 5, November 1998

F. Vető, Z. Horváth, T. Dóczi

**Biportal Endoscopic Management of Third Ventricle Tumors in
Patients with Occlusive Hydrocephalus: Technical Note**

Neurosurgery, Vol. 40, No. 4, 871-877, April 1997

J. A. Grotenhuis

Endoscope-Assisted Craniotomy

Techniques in Neurosurgery, Vol. 1, No. 3, 201-212, 1996

G. Fries, R. Reisch

**Biportal Neuroendoscopic Microsurgical Approaches to the
Subarachnoid Cisterns: A Cadaver Study**

Minimally Invasive Neurosurgery, Vol. 39, 99-104, 1996

A. Perneczky

Planning Strategies for the Suprasellar Region

Neurosurgeons 11, 343-348, 1992

AESCALAP NEUROSURGERY

Numerical Index

BN175R	89	FA066R	73	FD741B	52
		FA067R	61	FD743B	52
FA030R	74	FA067R	73	FD761B	52
FA031R	74	FA068R	61	FD763B	52
FA032R	74	FA068R	73	FD766B	53
FA033R	74	FA069R	85	FD767B	53
FA034R	74	FA070R	85	FD768B	53
FA035R	74	FA071R	85	FD769B	53
FA036R	74	FA072R	84	FD771B	50
FA037R	74	FA073R	84	FD771B	50
FA038R	74	FA074R	84	FD773B	50
FA039R	74	FA075R	84	FD774B	50
FA040R	74	FA076R	78	FD775B	50
FA041R	60	FD113D	62	FD776B	50
FA041R	72	FD114D	62	FD786B	53
FA042R	60	FD115D	62	FD787B	53
FA042R	72	FD116D	62	FD788B	53
FA043R	60	FD220R	83	FD789B	53
FA043R	72	FD222R	83	FD797B	54
FA044R	60	FD224R	83	FD798B	54
FA044R	72	FD226R	83	FD799B	54
FA045R	60	FD228R	82	FD805B	54
FA045R	72	FD228R	83	FD808B	54
FA046R	60	FD467R	55	FD809B	54
FA046R	72	FD701B	50	FD811B	54
FA047R	60	FD702B	50	FD812B	54
FA047R	72	FD703B	50	FD814B	54
FA060R	60	FD704B	50	FD815B	54
FA060R	72	FD705R	50	FD816B	54
FA061R	61	FD706B	50	FD818B	54
FA061R	73	FD711B	52	FD819B	54
FA062R	61	FD717B	51	FD821B	54
FA062R	73	FD718B	51	FD822B	54
FA063R	61	FD719B	51	FD823B	54
FA063R	73	FD731B	50	FD824B	55
FA064R	61	FD732B	50	FD825B	55
FA064R	73	FD733B	50	FD826B	55
FA065R	61	FD734B	50	FD827B	55
FA065R	73	FD735B	50	FD828B	55
FA066R	61	FD736B	50	FD831B	55

FD832B	55	FK906R	90	FM697R	44
FD833B	55	FK907B	90	FM700R	45
FD835B	55	FK907R	90	FM701R	45
FD836B	55	FK908B	90	FM702R	45
FD839B	55	FK908R	90	FM705R	45
FD840B	55	FK909B	90	FM706R	45
FD841B	55	FK909R	90	FM707R	45
FD848B	54	FK911B	64	FM710R	46
FD849B	54	FK912B	64	FM711R	46
FE489K	48	FK913B	64	FM712R	46
FE490K	48	FK914B	64	FM715R	46
FE491K	48	FK923B	90	FM716R	46
FE494K	48	FK923R	90	FM717R	46
FE495K	48	FK924B	90	FM720R	47
FE496K	48	FK924R	90	FM721R	47
FF151R	92	FK936R	90	FM722R	47
FF168R	92	FK937R	90	FM725R	47
FF280R	92	FK938R	90	FM726R	47
FF345R	80	FK966B	64	FM727R	47
FF357R	70	FK967B	64	FM730R	44
FF456B	63	FM156R	89	FM730R	45
FF457B	63	FM157R	89	FM730R	46
FF458B	63	FM158R	89	FM730R	47
FF459B	63	FM670R	44	FM731R	44
FF496R	65	FM671R	44	FM731R	45
FF496R	91	FM672R	44	FM731R	46
FF497R	65	FM675R	44	FM731R	47
FF497R	91	FM676R	44	FM732R	44
FF498R	65	FM677R	44	FM732R	45
FF498R	91	FM680R	45	FM732R	46
FF499R	65	FM681R	45	FM732R	47
FF499R	91	FM682R	45	FT489T	48
FH605SU	70	FM685R	45	FT490T	48
FH610R	71	FM686R	45	FT491T	48
FH611R	71	FM687R	45	FT494T	48
FH615	70	FM690R	44	FT495T	48
FK900B	64	FM691R	44	FT496T	48
FK901B	64	FM692R	44		
FK902B	64	FM695R	44	GA258SU	119
FK906B	90	FM696R	44	GA259SU	119

AESCULAP NEUROSURGERY

Numerical Index

GA261SU	119	GB698R	120	GF251R	87
GA262SU	119	GB941R	113	GF252R	87
GA263SU	119	GB942R	113	GF253R	87
GA264SU	119	GB943R	113	GF254R	87
GA702R	106	GB945R	113	GF255R	87
GA705R	106	GB947R	113	GF256R	87
GA706R	106	GF005B	58	GF257R	87
GA707R	106	GF025B	58	GF258R	87
GA708	106	GF026B	58	GF260R	87
GA710R	106	GF027B	58	GF261R	87
GA712R	106	GF030B	58	GF262R	87
GA713R	106	GF031B	58	GF263R	87
GA715R	106	GF032B	58	GF264R	87
GA722	116	GF033B	58	GF265R	87
GA724	116	GF035B	58	GF266R	87
GA736	117	GF036B	58	GF267R	87
GA737	118	GF037B	58	GF268R	87
GA749	113	GF038B	58	GF270R	87
GA761	107	GF039B	58	GF271R	87
GA762	107	GF045B	58	GF272R	87
GA763	107	GF046B	58	GF273R	87
GA764	107	GF047B	58	GF274R	87
GA800	104	GF048B	58	GF275R	87
GA806	104	GF049B	58	GF276R	87
GA808	104	GF055B	58	GF277R	87
GA822	116	GF056B	58	GF278R	87
GA824	116	GF057B	58	GF470R	59
GA836	117	GF058B	58	GF471R	59
GA837	118	GF059B	58	GF472R	59
GA849	113	GF235R	87	GF473R	59
GA861	105	GF240R	87	GF474R	59
GA862	105	GF241R	87	GF475R	59
GA863	105	GF242R	87	GF476R	59
GA864	105	GF243R	87	GF477R	59
GB600	120	GF244R	87	GF478R	59
GB600850	120	GF245R	87	GF479R	59
GB600860	120	GF246R	87	GF480R	59
GB691R	120	GF247R	87	GF481R	59
GB692R	120	GF248R	87	GK560R	88
GB697R	120	GF250R	87	GK580R	88

GK777R	57	GP135R	108	GP190R	112
GK780R	57	GP136R	108	GP193R	112
GK781R	57	GP141R	109	GP194R	112
GK785R	57	GP142R	109	GP195R	112
GK800R	88	GP143R	109	GP196R	112
GK801R	88	GP144R	109	GP201R	111
GK822R	56	GP145R	109	GP202R	111
GK823R	56	GP146R	109	GP203R	111
GK824R	56	GP147R	109	GP204R	111
GK825R	56	GP148R	109	GP205R	111
GK826R	56	GP149R	109	GP208R	111
GK827R	56	GP150R	109	GP209R	111
GK828R	56	GP151R	109	GP210R	111
GK829R	56	GP152R	109	GP210TC	112
GP111R	108	GP153R	109	GP211R	111
GP112R	108	GP154R	109	GP212R	111
GP113R	108	GP155R	109	GP213R	111
GP113TC	112	GP156R	109	GP214R	111
GP114R	108	GP157R	109	GP215R	111
GP115R	108	GP158R	109	GP216R	111
GP116R	108	GP161R	110	GP217R	111
GP117R	108	GP162R	110	GP218R	111
GP118R	108	GP163R	110	GP219R	111
GP119TC	112	GP164R	110	GP220R	111
GP120R	108	GP165R	110	GP224R	111
GP121R	108	GP166R	110	GP225R	111
GP122R	108	GP168R	110	GP226R	111
GP123R	108	GP169R	110	GP228R	111
GP124R	108	GP170R	110	GP229R	111
GP124TC	112	GP173R	110	GP230R	111
GP125R	108	GP174R	110	GP232R	112
GP126R	108	GP175R	110	GP233R	112
GP126TC	112	GP176R	110	GP234R	111
GP127R	108	GP177R	110	GP235R	111
GP128R	108	GP184R	108	GP236R	111
GP128TC	112	GP185R	108	GP238R	112
GP129R	108	GP186R	108	GP240R	112
GP129TC	112	GP187R	108	GP241R	112
GP133R	108	GP188R	112	GP242R	112
GP134R	108	GP189R	112	GP243R	112

AESCULAP NEUROSURGERY

Numerical Index

GP301R	114	GP353R	115	OK525R	80
GP302R	114	GP355R	115	OK560R	82
GP303R	114	GP356R	115	OK561R	82
GP304R	114	GP357R	115	OK562R	82
GP305R	114	GP358R	115	OK602R	78
GP306R	114	GP491R	117	OK603R	78
GP307R	114	GP492R	117	OK606R	78
GP308R	114	GP493R	117	OK607R	78
GP309R	114	GP494R	117	OK608R	78
GP311R	114	GP495R	117	OK609R	78
GP312R	114	GP496R	117	OK680R	79
GP313R	114	GP497R	117	OK681R	79
GP314R	114			OK682R	79
GP315R	114	JF324R	41	OK683R	79
GP316R	114	JG901	93	OK684R	79
GP317R	114	JG904	100	OP923	101
GP321R	114			OP940	100
GP322R	114	OF601R	89	OP940	101
GP323R	114	OK081R	77		
GP324R	114	OK082R	77	PE487A	71
GP325R	114	OK083R	77	PE506A	41
GP328R	114	OK084R	77	PE507A	71
GP329R	114	OK090R	76	PE526A	41
GP330R	114	OK090R	76	PV460	100
GP335R	115	OK091R	77	PV460	101
GP336R	115	OK092R	77	PV462	100
GP337R	115	OK093R	77	PV463	100
GP338R	115	OK094R	77	PV472	101
GP339R	115	OK105R	76	PV473	101
GP341R	115	OK106R	76	PV840	100
GP342R	115	OK107R	76	PV880	100
GP343R	115	OK108R	76	PV880	101
GP344R	115	OK505R	81	PV881	100
GP345R	115	OK506R	81	PV881	101
GP346R	115	OK507R	81	PV884	100
GP347R	115	OK508R	81	PV884	101
GP348R	115	OK509R	81	PV909	100
GP349R	115	OK520R	81	PV909	101
GP351R	115	OK521R	81	PV956	100
GP352R	115	OK522R	81	PV959	101

RT020R	93
RT040R	93
RT043R	93
RT044SU	93
RT046P	94
RT055P	94
RT060R	97
RT061R	97
RT063R	97
RT064R	97
RT065R	97
RT066R	97
RT068R	94
RT079205	94
RT079R	94
RT081R	94
RT090R	92
RT099R	70
RT099R	94
<hr/>	
TE730	119
TE734	119
TE734	119
TE734	119
TE780	119

AESFULAP[®] – a B. Braun brand

Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany

Phone +49 7461 95-0 | Fax +49 7461 95-2600 | www.aesculap.com

The main product trademark "Aesculap" is a registered trademark of Aesculap AG. Subject to technical changes. All rights reserved. This brochure may only be used for the exclusive purpose of obtaining information about our products. Reproduction in any form partial or otherwise is not permitted.