TECNOLOGIA  
PACCOCATH®

## AESCULAP® SeQuent® Please

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


Tecnologia SeQuent® Please

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
### First in man: ISR I/II study<sup>6,7</sup>

Treatment of coronary in-stent restenosis with a Paclitaxel coated balloon catheter - a prospective, randomized trial.

IN-STENT LATE LOSS AND RESTENOSIS AT 6 MONTHS AND MACE AT 2-YEARS:

	Uncoated balloon	Paclitaxel eluting balloon	p-value
In-stent late loss [mm]	0.81 ± 0.79	0.14 ± 0.46	0.001
In-segment late loss [mm]	0.80 ± 0.79	0.11 ± 0.48	0.001
In-stent restenosis	49%	6%	0.001
In-segment restenosis	51%	6%	0.001
MACE	46%	11%	0.001

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	SeQuent® Please	*PES (Paclitaxel eluting stent)
In-segment late loss [mm]	0.16 ± 0.38	0.49 ± 0.61
In-segment restenosis	5.5%	31.2%
MACE	6.1%	18.9%

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	SeQuent® Please	*PES	p-value
In-segment late loss [mm]	0.18 ± 0.41	0.39 ± 0.63	0.03
In-segment restenosis	6.7%	20.4%	0.03
MACE	4.7%	18.3%	<0.02

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## AESULAP® – a B. Braun brand

1- Axel DI et al. Paclitaxel Inhibits Arterial Smooth Muscle Cell Proliferation and Migration In Vitro and In Vivo Using Local Drug Delivery *Circulation* 1997;96:636-645 2- Hwang CW et al. Physiological Transport Forces Govern Drug Distribution for Stent-Based Delivery, *Circulation* 2001;104:600-605 3- Scheller B, Speck U et al. Paclitaxel Balloon Coating, a Novel Method for Prevention and Therapy of Restenosis *Circulation* 2004;110:810-814 4- Scheller B, Speck U, Böhm M Prevention of restenosis: is angioplasty the answer? *Heart* 2007;93:539-541 5- Speck U, Scheller B et al. Neointima Inhibition: Comparison of Effectiveness of Non-Stent Based Local Drug Delivery and a Drug eluting Stent in Porcine Coronary Arteries *Radiology* 2006;240:411-418 6- Scheller B et al. Treatment of Coronary In-Stent Restenosis with a Paclitaxel-Coated Balloon Catheter *New England Journal of Medicine* 2006; 355:2113-24 7- Scheller B, Hehrlein C, Bocksch W, Rutsch W, Haghi D, Dietz U, Böhm M, Speck U Two year follow-up after treatment of coronary in-stent restenosis with a paclitaxel-coated balloon catheter *Clin Res Cardiol* 97:773-781 (2008) 8- Drug-Eluting Balloons May Be Alternative to Drug-Eluting Stents *FDAnews Device Daily Bulletin* Oct. 26, 2007;Vol. 4;No. 211 9- Unverdorben M et al. in press Paclitaxel-Coated Balloon Catheter versus Paclitaxel-Coated Stent for the Treatment of Coronary In-stent Restenosis *Sequent® Please* - Registro nº 8.01369.90654

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B. Braun Brasil | S.A.C: 0800 0227286 | [www.bbraun.com.br](http://www.bbraun.com.br)