# IFU Lyoplast® Onlay

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# Lyoplant® Onlay

## Product Description

Lyoplant® Onlay is a lyophilized dura substitute of pure collagen derived from bovine pericardium and bovine split hide. The product is provided sterile, for single use only.

**Federal law restricts this device to sale by or on the order of a physician.**

## Composition

Lyoplant® Onlay consists of specially purified collagen. It is EtO sterilized.

## Indications for Use

Lyoplant® Onlay is indicated as a dura substitute for the repair of the dura mater.

## Contraindications

- Do not use in infected areas, open cranial trauma or open spina bifida.
- Do not use as a replacement for mechanically stressed connective tissue structures.
- Do not use as a substitute for parts in the arterial system or the cardiac wall.
- Do not use if there are known allergies to proteins of bovine origin.

## Precautions

Before using Lyoplant® Onlay, the surgeon should be familiar with the surgical technique, special applications and the properties of Lyoplant® Onlay in vivo.

## Warnings

- Lyoplant® Onlay should only be used if the package is intact. Visually inspect the packaging; do not use product if the package is found opened, punctured, torn or tampered with as sterility may be compromised.
- Do not re-sterilize, as the sterilization process may alter the structure of Lyoplant® Onlay and thus its properties in vivo.
- Do not use left-over implant pieces.
- Do not use product beyond expiration date.

## How Supplied

The product is supplied sterile and packaged as a single piece. The product is available in various dimensions ranging from 2.5 x 2.5 cm to 10.0 x 12.5 cm pieces. Specific sizes are listed by catalog number below:

Catalog Number	Dimensions
1067010	2.5 x 2.5 cm
1067020	5.0 x 5.0 cm
1067030	2.5 x 7.5 cm
1067040	7.5 x 7.5 cm
1067050	10.0 x 12.5 cm

## Possible Adverse Reactions

Hypersensitivity or an immune response to Lyoplant® Onlay may occur. The physician should closely monitor the patient for such a reaction and treat accordingly. Other possible adverse reactions may include, but are not limited to infection, postoperative adhesions, leakage of cerebrospinal fluid and hematoma.

## Recommendations for Use

To minimize tension after implantation, the correct size of Lyoplant® Onlay should be chosen and cut as closely as possible to fit the defect. Proper sizing of Lyoplant® Onlay is important for optimal results.

For improved flexibility and easier shaping, it is recommended that the graft be adequately rehydrated before use. Soak Lyoplant® Onlay in sterile 0.9% saline for a minimum of 3 minutes. Hydration is indicated by a color change from white (dry) to grey (wet).

Ensure that the fleece-like, porous surface is facing the dura (stamped "DURA SIDE" on graft). Identify the direction in which each surface should be facing before hydration.

Lyoplant® Onlay can be used as an onlay or sutured into place. When suturing, non-absorbable (polyester, polypropylene) suture material is recommended to improve suture security. The use of atraumatic round-bodied needles will allow suturing without significant damage to the implant.

## Storage

Store unopened product in a clean, dry and protected area between 20–30°C (68–86°F).

## Rev.

06/2013

# IFU Lyoplant® Onlay

## Intended use

Lyoplant® Onlay is used in neurosurgery as a dura substitute.

## Indications

Uses:

- For covering cerebral and cerebellar dura defects
- For decompressive craniectomies in elevated intracranial pressure
- For covering spinal dura defects
- For spinal decompression surgery

## Absolute contraindications

Do not use:

- In infected areas
- As a replacement for connective tissue structures that are subject to mechanical stress (tendons / ligaments)
- As a replacement for parts of the arterial system or the heart walls
- In cases of known foreign body sensitivity to implant materials; if fixation materials are used, the corresponding instructions for use must be followed
- For any application not listed under "Indications"

## Relative contraindications

Risks, side effects and interactions caused by the patient's comorbidities or in combination with other products are not known.

## Risks, side effects and interactions

The following risks, side effects and interactions are currently known:

- Infection
- Cerebrospinal fluid leakage
- Adhesion
- Foreign body reaction

Note:

The items in this list include the potential clinical consequences.

## Use

- Select a suitable implant size for the sealing of the defect.
- Cut Lyoplant® Onlay as required for the application situation. The edge of the implant should overlap the surrounding dura by approximately 1 cm.

**Note:** Rehydrate Lyoplant® Onlay in sterile physiological saline solution or another isotonic solution before implantation to improve smoothness and flexibility.

→ **Rehydration can be recognised by the colour change from white (dry) to grey (wet).**

**Note:** Lyoplant® Onlay consists of two different layers. It is important to ensure that the fleece-like porous side (labelled DURA SIDE) faces the dura. Check the orientation of the two layers before rehydration.

- Lay the implant flat on the defect margins, ensuring that it is not taut.
- Lyoplant® Onlay can be onlayed, or it can be sutured in place if required and considered necessary by the user (the use of non-absorbable suture material and atraumatic round-bodied needles is recommended).
- Fibrin glue can also be used for sealing.
- No other fixation methods may be used.

## Safety information and warnings

### a. User

- The following knowledge is required to successfully use this product:
  - Appropriate clinical training
  - Theoretical and practical command of all the surgical techniques required
  - In vivo properties of Lyoplant® Onlay
- Aesculap is not responsible for complications resulting from:
  - Incorrect indication, selection of incorrect implant size, or incorrect cutting of the implant
  - Using or combining the product other than as described in these instructions for use
  - Exceeding the limits of the treatment method or failure to observe basic medical precaution

### b. Product

- Safety with regard to the transmission of zoonosis. Based on the fact that bovine material from New Zealand is considered safe by European authorities with regard to BSE (bovine spongiform encephalopathy), the raw material is imported from this country. Additionally, Lyoplant® Onlay is treated with NaOH during further processing in order to further reduce the theoretical risk of transmission using this recognised depletion process.

### c. Information on shelf life, storage and reusability

- Damaged packages must not be used.
- Dispose of any implant pieces not required intraoperatively.
- If the implant is reused, the patient may suffer infection and the implant may lose its function. There are risks in the form of injury, illness or death from contamination and/or restricted functionality of the product.

**Danger**

**Danger to patients! Lyoplant® Onlay must not be resterilised because the structure of the implant and thus its in vivo behaviour may be adversely changed.**

► Danger due to changes in the material properties, Lyoplant® Onlay must not be resterilised!

- Lyoplant® Onlay must not be resterilised.
- Do not use Lyoplant® Onlay after the specified expiry date.
- Lyoplant® Onlay must be stored in a dry place at  $25 \pm 5$  °C.

## Product properties

Lyoplant® Onlay is a two-layered, three-dimensional implant made of purified collagen (type I/III) obtained from bovine pericardium and bovine split hide (composition:  $12 \pm 4$  mg/cm<sup>2</sup> pericardium;  $10 \pm 2$  mg/cm<sup>2</sup> spongy components). Lyoplant® Onlay is cleansed of non-collagenous components such as enzymes, lipids and non-collagenous proteins by the special preparation process. Lyoplant® Onlay is not chemically cross-linked.

The fleece-like properties of the porous side (bovine split hide) enables Lyoplant® Onlay to adhere to the defect margins.

Gentle freeze drying ensures that Lyoplant® Onlay retains its loose fibre structure, which promotes healing in the surrounding tissue after implantation.

After implantation, Lyoplant® Onlay is slowly degraded by enzymes and replaced by endogenous connective tissue. Connective tissue cells start to colonise the implant after just a few days. Complete revitalisation occurs within 1 to 3 months.

Lyoplant® Onlay has no independent pharmacological effect. The collagen implant is broken down into amino acids in the body.

## Delivery forms and package sizes

Lyoplant® Onlay is supplied EO sterilised for single use in double-packaged peel packs of various sizes.

## Single unit packaging

DIM	REF
2,5 x 2,5 cm	1067010
5,0 x 5,0 cm	1067020
2,5 x 7,5 cm	1067030
7,5 x 7,5 cm	1067040
10,0 x 12,5 cm	1067050

## Explanation of package symbols

	Follow the instructions for use
	Temperature limit
	Do not use if the package is damaged
	Attention
	Do not resterilise
	Keep dry
	Do not reuse
	MR safe
	Sterile as long as the package is undamaged and unopened. Sterilisation method: Ethylene oxide
	CE symbol and identification number of the listed location. The product meets the fundamental requirements of the council directive on medical products 93/42/EEC.
	Manufacturer
	Use by Year, Month, Day
	Date of Manufacture
	Batch identification
	Item number
	Size

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# Lyoplast® Onlay

## Description

Lyoplast® Onlay is an implant made from pure collagen obtained from bovine pericardium and bovine split-skin. In a special processing procedure, non-collagen components such as enzymes, fats and non-collagen proteins are removed from the Lyoplast® Onlay.

The gentle freeze-drying process ensures that the loose fiber architecture of Lyoplast® Onlay is preserved, which provides optimum integration and healing conditions after implantation.

After implantation, Lyoplast® Onlay is slowly broken down by enzymes and replaced by the body's own connective tissue. Lyoplast® Onlay has no pharmacological action of its own. The collagen implant is broken down by the body into amino acids.

## Composition

Lyoplast® Onlay is a biphasic three-dimensional collagen type I/III product. Lyoplast® Onlay is composed of  $12 \pm 4$  mg/cm<sup>2</sup> pericardium and  $10 \pm 2$  mg/cm<sup>2</sup> sponge components. Lyoplast® Onlay is not chemically cross-linked.

## Indication

Replacement and extension of connective tissue structures in neurosurgery:

- to close holes in the cerebral and cerebellar dura
  - after extirpation of tumors
  - to close CSF fistulas
- for procedures to relieve increased intracranial pressure
- to close holes in the spinal dura
  - after removal of spinal tumors
  - after spinal traumas
- for procedures to relieve spinal pressure

## Mode of application

The appropriate size of the Lyoplast® Onlay is chosen according to the area of application and cut according to the size of the defect, leaving an excess edge on the implant of approximately one centimeter (1 cm). Prior to implantation, Lyoplast® Onlay is placed in sterile saline solution or in another isotonic solution to ensure better suppleness and flexibility. The implant should be fixed without being pulled tautly. Lyoplast® Onlay consists of two different sides. Ensure that the fleece-like, porous side (labeled "DURA SIDE") is facing the dura. Identify the direction in which each side should be facing before hydration. Lyoplast® Onlay can be used as an Onlay graft. Ensure that the implant is lying flat against the defect edges and is not taut. In all other cases it is recommended to fixate the implant as well.

Lyoplast® Onlay can simply be applied, or, if required and if considered necessary by the user, it can also be sutured in place and sealed with fibrin glue.

Prior to implanting Lyoplast® Onlay, the user must be acquainted with the surgical technique, special applications and the characteristics of Lyoplast® Onlay in vivo.

## Mode of operation

After implantation, the Lyoplast® Onlay implant integrates with the body's own connective tissue cells, i.e. revitalization takes place. The fiber architecture of the compact side of the implant, the narrow thickness of the material, and a good connection with the recipient tissue are particularly important for this. The fleece-like quality of the porous side allows the Lyoplast® Onlay to adhere to the edges of the defect.

The implant starts to integrate with the connective tissue cells after just a few days. Full revitalization takes place within 1 to 3 months.

## Contraindications

Do not use in open cranial trauma or open spina bifida.

Lyoplast® Onlay should not be used

- in infected regions
- as a replacement for connective tissue structures that undergo mechanical stress (tendons/ligaments)
- as a replacement for parts of the arterial system or walls of the heart
- in case of known hypersensitivity to proteins of bovine origin.

## Warnings

**WARNING:** Risk of changing the material properties of Lyoplast® Onlay if resterilized!

Lyoplast® Onlay must not be resterilized as this can adversely alter the structure of the implant and therefore its behavior in vivo.

Lyoplast® Onlay may only be used if the packaging is undamaged.

Do not reuse already opened packages of Lyoplast® Onlay or implant pieces that were not required.

**WARNING:** Do not reuse the product!

Reuse of the implant can cause infections in the patient or user and functional loss of the implant. Risk of injury, illness or death through contamination and/or limited functionality of the product.

## Safety and precautions

Safety against the transmission of zoonoses.

Due to the fact that bovine material from New Zealand and Australia is considered BSE (bovine spongiform encephalopathy) -safe by the European authorities, the raw material is imported from these countries. Lyoplast® Onlay is also treated with sodium hydroxide solution during its further processing, and this recognized depletion process additionally reduces the theoretical risk of transmission.

## Side effects

- Occasionally, tissue adhesions can occur.
- Immune reactions cannot be ruled out with absolute certainty.
- The occurrence of CSF leaks cannot be fully ruled out in some individual cases.

## Sterilization

Lyoplast® Onlay is EO-sterilized.

## Other instructions

Do not use Lyoplast® Onlay after the expiry date indicated.

Use Lyoplast® Onlay immediately after opening the packaging.

Store Lyoplast® Onlay at  $25 \pm 5$  °C.

## Delivery forms and packaging sizes

Pack with 1 unit

DIM	REF
2,5 × 2,5 cm	1067010
5,0 × 5,0 cm	1067020
2,5 × 7,5 cm	1067030
7,5 × 7,5 cm	1067040
10,0 × 12,5 cm	1067050

## Revision

03/2016

## Explanation of the symbols on the packaging

	Follow the instructions for use
	Temperature limit
	Do not use if the package is damaged
	Attention
	Do not resterilise
	Keep dry
	Do not reuse
	MR-save
	Sterile as long as the package is undamaged and unopened. Sterilisation method: Ethylene oxide
	CE symbol and identification number of the listed location. The product meets the fundamental requirements of the council directive on medical products 93/42/EEC.
	Manufacturer
	Use by Year, Month, Day
	Date of Manufacture
	Batch identification
	Item number
	Size

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Aesculap AG | Am Aesculap-Platz | 78532 Tuttlingen | Germany