

Aesculap[®] Spine Classics

Micro Lumbar Discectomy System



Aesculap Spine

Spine Classics

■ Minimally invasive approach to the posterior lumbar spine

■ Optimized basic set

composed by experts with a longterm experience in thousands of surgeries

■ Autoclavable tray

efficient storage of the instruments

■ Rapid ball snap closure

easy and fast handling secure connection compatibility of all CASPAR blades and retractors

■ Non-slip-fenestrated retractors

secure fit in muscles

■ Special retractors and blades for the treatment of extraforaminal disc herniations

'classic' medial and 'special' lateral paraspinal approach

■ Indications


- Lumbar fusion
- PLIF
- TLIF
- Micro discectomy
- Laminectomy, Foraminotomy



Example of use



A) Lumbar minimally invasive approach for discectomy



B) Minimally invasive approach to the lumbar spine for a posterior lumbar fusion



Spine Classics

Blades, medial,
surface black-coated

Ref. No.	Length of Blade
BV851R	40 mm
BV852R	45 mm
BV853R	50 mm
BV854R	55 mm
BV855R	60 mm
BV856R	65 mm
BV857R	70 mm



Blades, lateral, with ball snap closure,
surface black-coated

Ref. No.	Length of Blade
BV870R	35 mm
BV871R	40 mm
BV872R	45 mm
BV873R	50 mm
BV874R	55 mm
BV875R	60 mm
BV876R	65 mm
BV877R	75 mm
BV878R	85 mm



**Specula,
surface black-coated**

Ref. No.	Length of Blade
BV799R	40 mm
BV800R	45 mm
BV801R	50 mm
BV802R	55 mm
BV803R	60 mm
BV804R	65 mm
BV806R	75 mm
BV808R	85 mm



**Counter-Retractor, with ball snap closure,
surface black-coated**

Ref. No.
BV812R



**Forceps for changing blades with
ball snap closure**


Ref. No.
BV399R



Spine Classics

Spine Classics MLD Basic-Instruments

Ref. No.	Description
BV071R	Self-Retaining Retractor, Length of Blades 13 mm, sharp
DO439R	Retractor, 55 x 14 mm
BT465R	Retractor, 67 x 17 mm
FF699R	Nerve Root Retractor, 6 mm wide
FD347R	Hook Knife
FF307R	Micro-Dissector, 2 mm, downwards curved
BT090R	Exploration and Coagulation Hook, probe end, non insulated
FK703R	Scoop, straight, 4.4 mm
FK713R	Scoop, angled, 4.4 mm
FF532R	Rongeur, straight, 160 x 2.0 mm
FF533R	Rongeur, straight, 160 x 3.0 mm
FF534R	Rongeur, straight, 160 x 4.0 mm
FF543R	Rongeur, up bite, 160 x 3.0 mm
FF544R	Rongeur, up bite, 160 x 4.0 mm
FF523R	Rongeur, down bite, 160 x 3.0 mm
FF717R	Bone Punch, upwards cutting, 2.0 mm
FF720R	Bone Punch, upwards cutting, 3.0 mm
FF724R	Bone Punch, upwards cutting, 4.0 mm
FF728R	Bone Punch, upwards cutting, 5.0 mm with over-size jaw opening
JF223R	Wire Basket 540 x 253 x 64 mm
JF932	Silicon Mat
JF120R	Storage Frame 240 x 225 x 161 mm



Spine Classics

FG170 Spine Classics MLD Retractor Set

Ref. No.	Description
BV799R	Speculum 40 mm
BV800R	Speculum 45 mm
BV801R	Speculum 50 mm
BV802R	Speculum 55 mm
BV803R	Speculum 60 mm
BV804R	Speculum 65 mm
BV812R	Counter Retractor, only
BV870R	Blade 35 mm
BV871R	Blade 40 mm
BV872R	Blade 45 mm
BV873R	Blade 50 mm
BV874R	Blade 55 mm
BV875R	Blade 60 mm
BV876R	Blade 65 mm
BV399R	Forceps for changing Blades with Ball Snap Closure
FG171P	Spine Classics MLD Tray I, only
FG172P	Spine Classics MLD Tray II, only
JF223R	Wire Basket 540 x 253 x 64 mm



